

Product Information

Anti-Luciferase

produced in rabbit, IgG fraction of antiserum

Catalog Number **L0159**

Product Description

Anti-Luciferase is produced in rabbit using as immunogen firefly (*Photinus pyralis*) luciferase. Whole antiserum is purified to provide an IgG fraction of antiserum.

Analysis of gene expression is most commonly assayed by transient transfection. These systems are usually based on the use of fusion genes which are inserted into cells, and expression of the gene is assayed within 48 hours after introduction of DNA. Usually the fusion consists of the promoter binding site or enhancer sequence under study which is attached to a reporter gene. The amount of the reporter protein synthesized under the experimental conditions is presumed to reflect the ability of the sequences studied to direct or promote transcription. Several enzymes are commonly used as reporter proteins, among them are chloramphenicol acetyl transferase (CAT), β -galactosidase, human growth hormone (hGH), and luciferase. Luciferase has become one of the more widely used reporter enzymes. The reporter plasmid contains the gene from the firefly *Photinus pyralis*. The enzyme catalyzes a bioluminescent reaction which requires the substrate luciferin as well as Mg^{+2} and ATP. Mixing these reagents with the cell extract containing luciferase results in a flash of light that decays rapidly. This light can be detected by a luminometer. The total light emission is proportional to the luciferase activity of the sample. The luciferase assay is fast and sensitive and does not require a radioactive substrate as in the CAT assay. A disadvantage of the luciferase assay is that it requires a rather expensive instrument, the luminometer, to measure enzyme activity. In addition, this assay lacks reproducibility between samples, largely due to the rapid kinetics of the emission.

Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, with 15 mM sodium azide as a preservative.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for a maximum of one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Product Profile

The use of an antibody to detect luciferase can provide an alternative detection assay. This assay directly detects luciferase protein levels, and thus has the advantage that it does not require luciferase activity and is not dependent on rapid kinetics. Moreover, antibodies can detect the luciferase enzyme expression *in situ*, providing a means to study the localized signal sequences using luciferase as a reporter gene.

Immunofluorescence: a working dilution of 1:500-1:1,000 is obtained on methanol-acetone fixed transfected* cells using an FITC conjugated secondary antibody.

* Transfected with a reporter plasmid containing the gene for luciferase.

Immunoblotting: a working dilution of 1:10,000 is determined using purified luciferase

Note: To obtain best results, it is recommended that each user determine the optimal working dilution for individual applications by titration assay.

References

1. De Wet, J. R., Wood, K.V., et al., *Mol. Cell Biol.* **7**, 725 (1987).

DS,PHC 08/15-1