1.14687.0001

Spectroquant® Crack Set 10

1. General

When performing photometric measurements, as a rule only free ions are determined. In wastewater and in other matrices, however, above all heavy-metal ions are present in complex-bound form. As a measure to determine the **total content** of dissolved metal - i. e. also the complex-bound portion - it is thus necessary to perform a digestion when preparing the sample. This is the purpose of the Crack Set 10 in connection with a thermoreactor.

The Crack Set 10 can also be used for sample pretreatment for the determination of the total phosphorus content.

2. Method

Heavy-metal complexes and phosphorus-containing compounds are destroyed by digestion with sulfuric acid and peroxodisulfate. Chromium(III) is oxidized to form chromium(VI), phosphorus compounds are converted into orthophosphate.

3. Number of digestions

The reagents are sufficient to conduct 100 digestions.

4. Applications

Sample pretreatment for the determination of cadmium, chromium, copper, iron, lead, nickel, total phosphorus, and zinc

**The COD of the sample may not exceed 200 mg/l.
The addition of acid to cyanide-containing samples can result in the formation of toxic hydrocyanic acid!**

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

1 bottle of reagent R-1
1 bottle of reagent R-2
1 bottle of reagent R-3
1 pack of universal indicator paper pH 1 - 14 (Cat. No. 1.10962)
1 blue dose-metering cap

Accessories:

Empty cells 16 mm with screw caps (25 pcs), Cat. No. 1.14724
Thermoreactor
Pipette for a pipetting volume of 10 ml

6. Preparation

At the first use **replace the screw cap of the reagent bottle R-2 by the blue dose-metering cap.**

Hold the reagent bottle **vertically** and, at each dosage, press the slide **all the way** into the dose-metering cap. **Before each dosage** ensure that the slide is **completely retracted.**

Reclose the reagent bottle with the screw cap at the end of the measurement series, since the function of the reagent is impaired by the absorption of atmospheric moisture.

7. Procedure

Pretreated sample (pH 5 - 8)	10 ml	Pipette into an empty cell.
Reagent R-1	1 drop ¹⁾	Add and mix.
Reagent R-2	1 dose	Add, close the cell tightly , and mix.

Heat the cell at 120 °C in the preheated thermoreactor for 1 hour. Allow the closed cell to cool to room temperature in a test-tube rack.

Do not cool with cold water!

Reagent R-3 ²⁾	3 drops ¹⁾	Add to the cool cell and mix.
---------------------------	-----------------------	-------------------------------

Check the pH of the solution with universal indicator paper. Adjust the pH, if necessary, according to the instructions given in the package insert of the test kit in question.

¹⁾ **Hold the bottle vertically while adding the reagent!**

²⁾ In the case that the iron concentration is to be measured after the digestion, **the addition of R-3 is not necessary.**

The digestion solution can be analyzed directly with the test kits given in the table:

Parameter	Cat. No.	Parameter	Cat. No.
Cadmium	1.14834 ²⁾	Lead	1.14833 ²⁾
	1.01745		1.09717
Chromate	1.14758	Nickel	1.14554 ²⁾
Cobalt	1.17244 ²⁾		1.14785
Copper	1.14553 ²⁾	Phosphate	1.14848
	1.14767	Zinc	1.14566 ²⁾
Iron ¹⁾	1.14549 ²⁾		1.00861 ²⁾
	1.14896 ²⁾		1.14832
	1.00796		
	1.14761		

¹⁾ **When analyzing iron, it is not necessary to add reagent R-3!**

²⁾ Cell tests

8. Notes

- Reclose the reagent bottles immediately after use.
- Information on disposal can be obtained at www.disposal-test-kits.com.**

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