

## Product Information

### **Pepsin** from porcine gastric mucosa

Product Number **P 6887**  
Storage Temperature -0 °C

#### **Product Description**

MW: 34,620 (amino acid sequence)<sup>1</sup>  
pI: 2.2 - 3.0<sup>2</sup>; 2.2, 2.8<sup>3</sup>  
 $\lambda_{\max}$ : 278 nm<sup>4</sup>  
Extinction coefficient:  $E^M = 51,300$ <sup>4</sup>

Pepsin, unlike some other peptidases, hydrolyzes only peptide bonds, not amide or ester linkages. The cleavage specificity includes peptides with an aromatic acid on either side of the peptide bond, especially if the other residue is also aromatic or a dicarboxylic amino acid. Increased susceptibility to hydrolysis occurs if there is a sulfur-containing amino acid close to the peptide bond, which has an aromatic amino acid. Pepsin will also preferentially cleave at the carboxyl side of phenylalanine and leucine and to a lesser extent at the carboxyl side of glutamic acid residues. Pepsin will not cleave at valine, alanine, or glycine linkages.<sup>5</sup> Some good substrates of pepsin are Z-L-tyrosyl-L-phenylalanine, Z-L-glutamyl-L-tyrosine, or Z-L-methionyl-L-tyrosine. Amidation of the C-terminal carboxyl group prevents hydrolysis by pepsin.<sup>5,6</sup>

Pepsin is commonly used in the preparation of Fab fragments from antibodies. The optimal pH for the pepsin reaction is 1.5-2.5, which will not be detrimental to the antibody, if it is not exposed for long durations to the low pH. Solutions should be adjusted to neutral pH for storage. The control of pepsin digestion of antibodies has been reported.<sup>7</sup>

For general digestion of proteins, suggested conditions are a 0.4% solution of pepsin in 10 mM HCl, and digestion for 30-90 minutes at 37 °C. Pepsin has optimal activity with native proteins at approximately pH 1.0, but with some denatured proteins the optimal activity is at approximately pH 1.5-3.5.<sup>8,9</sup>

Pepsin is inhibited by several phenylalanine-containing peptides.<sup>10</sup>

#### **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

#### **Preparation Instructions**

Pepsin is soluble in deionized water at 1% (10 mg/ml) and at 0.4% (4 mg/ml) in cold 10 mM hydrochloric acid.

#### **Storage/Stability**

Solutions at pH 4.4 are stable at -20 °C for about 2-3 months.<sup>11</sup> Pepsin is not active when not at an acidic pH and a solution is stable at pH 6-7. Bringing the pH up to 8; however, will irreversibly inactivate pepsin. Pepsin is irreversibly denatured at pH 8.5 - 11 at room temperature.<sup>12</sup>

#### **References**

1. Sepulveda, P., et al., Primary Structure of Porcine Pepsin. III. Amino Acid Sequence of a Cyanogen Bromide Fragment, CB2A, and the Complete Structure of Porcine Pepsin. *J. Biol. Chem.*, **250**, 5082 (1975).
2. Jonsson, M., Isoelectric Spectra of Native and Base Denatured Crystallized Swine Pepsin. *Acta Chem. Scand.*, **26**, 3435-3440 (1972).
3. Malamud, D., and Drysdale, J.W., Isoelectric Points of Proteins: A Table, *Anal. Biochem.*, **86**, 620-647 (1978).
4. *Proc. Natl. Acad. Sci.*, **45**, 915-922 (1959).
5. Sweeney, P.J., and Walker, J.M., in *Enzymes of Molecular Biology*, Burrell, M.M., ed., Humana Press (Totowa, NJ: 1993), pp. 290-291.
6. *Enzymes*, Dixon, M., et al., Academic Press (New York, NY: 1979), p. 262.
7. Rea, D.W., and Ultee, M.E., A Novel Method for Controlling the Pepsin Digestion of Antibodies. *J. Immunol. Meth.*, **157**, 165-173 (1993).
8. *Arch. Biochem. Biophys.*, **57**, 163-173 (1955).
9. *J. Biol. Chem.*, **234**, 3137-3145 (1959).

10. Knowles, J.R., et al., The pH-dependence of the Binding of Competitive Inhibitors to Pepsin. *Biochem. J.*, **113**, 343-51 (1969).
11. Rajagopalan, T.G., et al., Pepsin from Pepsinogen. Preparation and Properties. *J. Biol. Chem.*, **241**, 4940 (1966).
12. Ryle, A.P., The Porcine Pepsins and Pepsinogens. *Methods in Enzymol.*, **19**, 316-336 (1970).

ARO 9/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.