



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

(+,-)-THIOPENTAL SODIUM/SODIUM CARBONATE DEA SCHEDULE III Product Number **T6023**

CAS #: 71-73-8

Synonyms for Thiopental (free acid): 5-Ethyl-5-(1-methylbutyl)-2-thiobarbituric acid; Pentothal; Penthiobarbital; Thionembutal; Thiopentone

Product Description

Appearance: Light yellow powder with a garlic-like odor

Molecular formula: $C_{11}H_{17}N_2O_2SNa$

Formula weight: 264.3

Dissociation Constant: $pK_a = 7.6$ at $20^\circ C$

The name thiopental sodium has been used to describe preparations of both thiopental sodium alone or physical mixtures of thiopental sodium and sodium carbonate (buffer salt).¹ This product (T6023) is a mixture of approximately 100 parts sodium thiopental and 6 parts sodium carbonate. Sodium thiopental is a barbiturate used as a short-term intravenous anesthetic. It was first introduced for surgical use in 1935 by Lundy.⁶ The anesthetic effect is due to initial distribution of the material to the brain. Thiopental will then redistribute to other tissues, primarily body fat.^{1,4,6} The barbiturate significantly binds to plasma proteins and it is metabolized by the liver, where the major pathway appears to be side-chain oxidation.^{4,6} Anesthetic IV doses for animal research fall in the range of 20-30 mg/kg depending on the species.⁷ Several references to analytical procedures have been published.^{2,4,5}

Preparation Instructions

Sigma tests the solubility of this mixture at 50 mg/ml in water and obtains a clear, faint yellow to yellow-green solution. Solutions are alkaline (a pH range of 10.2-11.2 for a 8% solution) and decompose upon standing.³ A 2.5% solution in water was stable for at least 10 days at $25^\circ C$.⁴ Upon boiling a precipitate will form.² Dispose of solutions if they become cloudy or a precipitate forms.⁴ In an acidic solution, the material will convert to the free acid which is water insoluble and will precipitate.^{1,4} Buffers containing Ca^{2+} or Mg^{2+} should be avoided due to formation of insoluble carbonate salts, as well as work with amine salts that liberate the free base in alkaline solutions.⁴

Storage/Stability

The powdered product is hygroscopic. It should be stored desiccated at room temperature.

References

1. Martindale: The Extra Pharmacopoeia, 31th ed., p. 1265.
2. United States Pharmacopeia.
3. Merck Index, 12th ed., #9487 (1996).
4. McLeish, M.J., Anal. Prof. Drug. Subst. Excip., 21, 535-572 (1992).
5. Clarkes Isolation and Identification of Drugs, p. 1017.
6. Goodman and Gilman's, The Pharmacological Basis of Therapeutics, Eighth edition, p. 301.
7. Drug Dosage in Laboratory Animals, A Handbook, 3rd Ed., p. 471.

feb 4/27/00

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.