



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sia.com  
sigma-aldrich.com

## Product Information

### Glycerol Cell Culture Tested

Product Number **G 2025**  
Store at Room Temperature

#### Product Description

Molecular Formula:  $C_3H_8O_3$   
Molecular Weight: 92.1  
CAS Number: 56-81-5  
Melting point: 17.8 °C  
Density: 1.26 g/ml  
Synonyms: glycerin, 1,2,3-propanetriol

This product is cell culture tested and insect cell culture tested. It is appropriate for use in cell culture and insect cell culture experiments.

Glycerol is a viscous liquid that is widely used in molecular biology research. Applications include purification of bacteriophage  $\lambda$  particles by step gradients, calcium-phosphate-mediated transfection of eukaryotic cells with plasmid DNA, DEAE-dextran-mediated transfection of DNA into cells, DNA sequencing reactions, and mediation of mispriming in PCR. Glycerol is used in the preparation of DNA and protein samples for both agarose and polyacrylamide gel electrophoresis. At 5-10% concentration, glycerol increases the density of a sample so that the sample will layer at the bottom of a sample well on the electrophoretic gel. Glycerol is also used in casting gradient gels.<sup>2</sup>

The use of glycerol in improved focusing of alkaline pH protein samples in immobilized pH gradient (IPG) strips has been published.<sup>3</sup> Studies of the effect of glycerol on separation of DNA fragments by capillary electrophoresis have been reported.<sup>4</sup>

Glycerol is used in the long-term stabilization of enzyme solutions and storage of bacterial stocks.<sup>2</sup> The specific gravity and freezing point values for various aqueous solutions of glycerol have been reported.<sup>1</sup>

#### Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

#### Preparation Instructions

Glycerol is miscible with water and alcohol at all concentrations<sup>1</sup>.

#### References

1. The Merck Index, 12th ed., Entry# 4493.
2. Molecular Cloning: A Laboratory Manual, 3rd ed., Sambrook, J. F., et al., Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY: 2001), pp. 2.52-2.53, 8.9, 8.23, 8.78, 12, 108-12.109, 13.90, 16.14-16.15, 16.28.
3. Hoving, S., et al., Preparative two-dimensional gel electrophoresis at alkaline pH using narrow range immobilized pH gradients. *Proteomics*, **2(2)**, 127-134 (2002).
4. Liang, D., et al., Effect of glycerol-induced DNA conformational change on the separation of DNA fragments by capillary electrophoresis. *J. Chromatogr. A*, **931(1-2)**, 163-173 (2001).

GCY/NSB 2/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.