

Analytix *Notes*

 **Fluka**
Riedel-deHaën®

Disks and Strips for Microbiology



SIGMA-ALDRICH

Disks and Strips for Microbiology Disks and Strips are very helpful for identification and confirmation of microorganisms or to monitor sterilization. They are based on rapid methods, are easy to prepare and have a very good price/quality relationship.

The strips and disks are made of a cellulose based material impregnated with the appropriate chemical reagents. This converts them into intelligent systems that can be used for the detection of specific abilities and properties of microorganisms, either based on the detection of enzymes using chromogenic substrates or on complex building reactions. Also, sensitivity to certain inhibitory substances can be also be tested for.

Contact information

Jvo Siegrist
Product Manager Microbiology
Fluka Riedel de Haën
Phone: ++41817552449
Fax: ++41817552824
E-mail: isiegris@europe.sial.com
www.sigma-aldrich.com/microbiology

- | | |
|---|--|
| <p>03 <i>Aminopeptidase Test</i>
(Gram-Positive-Test)
Cat. No. 75554</p> <p>04 <i>Bacitracin Disks</i>
Cat. No. 08382</p> <p>04 <i>Bile Esculin Disks</i>
(Esculin Bile Disks)
Cat. No. 80507</p> <p>05 <i>DMACA Indole Disks</i>
(Indol Detection Disks)
Cat. No. 05686</p> <p>06 <i>Hippurate Disks</i>
Cat. No. 40405</p> <p>07 <i>Hydrogen Sulfide Test Strips</i>
(Lead acetate test strips, H₂S test strips)
Cat. No. 06728</p> <p>08 <i>Indoxyl Strips</i>
(Acetoxyindol Strips)
Cat. No. 04739</p> <p>09 <i>Kovac's Reagent Strips</i>
(Indole Reagent Test strips according to Kovac)
Cat. No. 78719</p> <p>10 <i>β-Lactamase Strips</i>
Cat. No. 56348</p> <p>11 <i>Nitrate Reagent Disks</i>
Cat. No. 08086</p> <p>12 <i>Nitrocefin Disks</i>
Cat. No. 49862</p> | <p>13 <i>ONPG Disks</i>
(2-Nitrophenyl β-D-galactopyranoside Disks,
β-Galactosidase Test Disks)
Cat. No. 49940</p> <p>14 <i>Optochin Disks</i>
Cat. No. 74042</p> <p>14 <i>Oxidase Strips</i>
Cat. No. 40560</p> <p>16 <i>PYRase Strips</i>
(Pyrrolidonyl Peptidase Strips)
Cat. No. 67886</p> <p>17 <i>Oxidase Test Disks</i>
Cat. No. 70439</p> <p>18 <i>Sterility Indicator</i>
(Steam Sterilization)
Cat. No. 74041</p> <p>18 <i>Sterile Disks</i>
Cat. No. 74146</p> <p>19 <i>Sterility Indicator</i>
(Radiation Sterilization)
Cat. No. 05290</p> <p>20 <i>Tributyrin-Strips</i>
(TRIBU Strips)
Cat. No. 75744</p> <p>21 <i>X Factor Differentiation Disks</i>
Cat. No. 77148</p> <p>21 <i>V Factor Differentiation Disks</i>
Cat. No. 89788</p> <p>21 <i>X + V Factors Differentiation Disks</i>
Cat. No. 8482</p> <p>22 <i>Antimicrobial Susceptibility Disks</i></p> |
|---|--|

Fluka Cat. No. 75554 Aminopeptidase Test (Gram-Positive Test)

The aminopeptidase test is intended for use in the detection of Gram-positive microorganisms by checking for the presence of L-alanine aminopeptidase.

Composition

Each package contains 50 test strips.
The reaction zone of each strip contains 0.5 µmol L-alanine-4-nitroanilide and buffering agents.

Storage

Store in a dry place.

Instructions

Remove a small sample from a single colony with an inoculation loop and suspend it in 0.2 mL distilled water placed in a test tube. Place a test strip into the clear opalescent suspension and incubate for between 10-30 minutes at 37°C. In presence of L-aminopeptidase, the solution becomes yellow, indicating that the microorganisms are

Gram-negative. (*Bacteroides vulgatus*, *Bacteroides fragilis*, *Campylobacter sp.*, *Veillonella parvula* are exceptions to this general rule.) If no yellow colouration appears, then L-alanine aminopeptidase is absent and therefore the microorganisms are Gram-positive.

Note: The growth medium from which the colonies were taken should not contain any dyes or indicators. The use of pigmented colonies is not recommended .

Principle and Interpretation of Results

L-alanine aminopeptidase is an enzyme from the bacterial cell wall that cleaves L-alanine from various peptides and it is found almost only in Gram-negative microorganisms. Gram-positive or Gram-variable microorganisms show no or very weak activity. The aminopeptidase test is a reliable method for determining Gram behaviour, however, it does not replace Gram-staining, as it cannot show morphology.

Organisms

Gram-negative bacteria
(Exceptions: *Bacteroides vulgatus*, *Bacteroides fragilis*, *Campylobacter sp.*, *Veillonella parvula*)
Gram-positive bacteria

L-alanine aminopeptidase

Present (yellow)
Not present (no coloration)

References


- (1) G.M. Carlone, M.J. Valdez, and M.J. Pickett. Method for Distinguishing Gram-positive from Gram-negative Bacteria., *J. Clinical Microbiology*, 16,1157 (1982)
- (2) E.H. Lennette, A. Balows, W.J. Haulser, and J.P. Truant (eds.), *Manual of Clinical Microbiology*, 3rd edition. American Society for Microbiology. (1980)
- (3) J.D. Costin, M. Kappner, W. Schmidt: Differenzierung von Gram-positiven Bakterien und Gram-negativen Bakterien mit dem L-Alanin Aminopeptidase Test, *Forum Mikrobiol.*, 351 (1983)
- (4) G. Cerny, Method for Distinction of the Gram-Negative from Gram-Positive Bacteria, *Eur. J. Appl. Microbiol.*, 3, 223 (1976)
- (5) G. Cerny, Studies on the Aminopeptidase-Test for the Distinction of the Gram-Negative from Gram-Positive Bacteria, *Eur. J. Appl. Microbiol. Biotechnol.*, 5, 113 (1978)
- (6) I. Otte, A. Tolle, Aminopetidase- und Gram-Reaktion von Bakterien, *Milchwiss.*, 35, 215 (1980)

Picture

Aminopeptidase Test

(from left to right)

1. Negative Reaction
2. Postive Reaction
3. Postive Reaction



Fluka Cat. No. 08382 Bacitracin Disks The bacitracin disks are used in the presumptive identification of group A β -hemolytic *Streptococci* and allow the differentiation of group A β -hemolytic *Streptococci* from other β -hemolytic *Streptococci*.

The bacitracin disks are used in the presumptive identification of group A β -hemolytic *Streptococci* and allow the differentiation of group A β -hemolytic *Streptococci* from other β -hemolytic *Streptococci*. The bacitracin test should be performed together with the SXT susceptibility test (Fluka Cat. No. 73477), as their combined results increase the sensitivity and the predictive value of the bacitracin test.

Composition

Each package contains 50, 6 mm diameter sterile filter paper disks impregnated with 0.04 Units of bacitracin.

Storage

Store in the freezer below -18°C in the containers provided. Allow to equilibrate to room temperature before opening. Return to freezer storage immediately after use.

Instructions

Prepare blood agar (Fluka Cat. No. 70133) plates and inoculate the plates with the suspect organism using surface spreading technique to obtain confluent growth. Place the bacitracin disks aseptically onto the inoculated surface and press gently. Invert the plates and incubate at $35-37^{\circ}\text{C}$ in 5-10 % CO_2 for 18-24 hours, until colony growth is observed. Check for the zone of inhibition around the disk.

Interpretation of results

A zone of inhibition greater than or equal to 14 mm indicates susceptibility to bacitracin and is presumptive of group A *Streptococci*. For further confirmation serological grouping is recommended.

Quality control


The table below illustrates a range of performance control strains in routine use (blood agar plate incubated at $35-37^{\circ}\text{C}$ for 18-24 hours).

Test Organisms (ATCC)	Diameter of zone of inhibition
<i>Streptococcus pyrogenes</i> (19615)	15 mm
<i>Streptococcus agalactiae</i> (27956)	< 14 mm

References

- (1) E.M. Barnes, G.C. Mead, C.S. Impey, B.W. Adams, The Effect of Dietary Bacitracin on the Incidence of *Streptococcus faecalis* Subspecies Liquefaciens and Related *Streptococci* in the Intestines of Young Chicks. Brit Poult Sci 19: 713-723. (1978)
- (2) E.J. Baron, S.M. Finegold, Bailey and Scott's Diagnostic Microbiology, 8th ed. St. Louis: Mosby (1990)
- (3) A. Balows, W.J. Hausler, K.L. Herman, et al., Manual of Clinical microbiology, 5th ed. Washington, DC: ASM (1991)
- (4) MacFaddin J.F., Biochemical Tests for the Identification of Medical Bacteria., 3rd ed. Philadelphia: Lippincott Williams & Wilkins (2000)

Fluka Cat. No. 80507 Bile Esculin Disks (Esculin Bile Disks) The bile esculin disks are used for rapid detection of esculin hydrolysis in presence of bile.



Picture Bile Esculin Disks

The detection of esculin hydrolysis in presence of biles allows the differentiation of group D *Streptococci* from non-group D *Streptococci*, as only group D *Streptococci* can hydrolyze the esculin to esculentin and dextrose. The resulting esculentin then reacts with iron salts such as ferric citrate, being a blackish-brown coloured complex formed.

Composition

Each package contains 50, 6mm diameter sterile filter paper disks impregnated with esculin.

Instructions

Place a bile esculin disk on the seeded bile esculin agar base (without esculin) plate or another media. Incubate at 35°C for 18-24 hours.

Quality control

Culture characteristics after 18-24 hours at 35°C.

References

- (1) Rochaix, C.R.Soc. Biol., 90, 771 (1924)
- (2) Meyer and Schönfeld, Zentralbl. Bacteriol. Parasitenkd. Infektionskr. Hyg. Abt. I Orig., 99, 402 (1924)

- (3) J.F. MacFaddin, Biochemical Tests for Identification of Medical Bacteria, 2nd ed., Williams and Wilkins, Baltimore (1980)
- (4) A.E. Greenberg, R. R. Trussell and L. S. Clesceri (Eds.), Standard Methods for the Examination of Water and Wastewater, 16th ed., A.P.H.A., Washington D.C. (1985)
- (5) R.R.Facklam, M.D. Moody, Presumptive Identification of Group D *Streptococci*: the Bile-Esculin Test. Appl. Microbiol., 20, 245 (1970)
- (6) S.C. Edberg, S. Pittman, J.M. Singer, Esculin Hydrolysis by *Enterobacteriaceae*., J. Clin. Micro. 6, 111 (1977)

Test Organisms (ATCC)

<i>Streptococcus faecalis</i> (29212)
<i>Streptococcus pyrogenes</i> (19615)
<i>Listeria monocytogenes</i> (19118)

Esculin Hydrolysis

+
-
+

Fluka Cat. No. 05686 DMACA Indole Disks (Indol Detection Disks)

The DMACA indole disks are used to determine the ability of an organism to split tryptophan into indole and α -aminopropionic acid.

The DMACA indole disks are used to determine the ability of an organism to split tryptophan into indole and α -aminopropionic acid. The presence of indole can be detected by adding DMACA, that gives rise to a blueish-purple complex. Using this method it is possible to differentiate *Escherichia coli* from *Klebsiella*.

Contents

Each package contains 50 disks impregnated with p-dimethylaminocinnamaldehyde (DMACA).

Quality control

Test Organisms (ATCC)

<i>Escherichia coli</i> (25922)
<i>Pseudomonas aeruginosa</i> (27853)
<i>Klebsiella pneumoniae</i> (13883)

Instructions

Place the disk on a suspect colony grown in a media such as modified HiCrome UTI agar (Fluka Cat. No. 16636) or Christensen's urea agar (Fluka Cat. No. 27048). Lay the Petri dish against a white background. Observe for the appearance of blue-purple colour within 10-30 seconds.

DMACA

+
-
-

References

- (1) R. Vracko, J.C. Sherris, Indole-spot test in bacteriology., Am. J. Clin. Pathol., 39, 429 (1963)
- (2) V.L. Sutter, W.T. Carter, Evaluation of Media and Reagents for Indole-Spot Test in Anaerobic Bacteriology., Am. J. Clin. Pathol., 58, 335 (1972)
- (3) G.D. Fay, A.L. Barry, Methods for Detecting Indole Production by Gram-negative Non-Spore Forming Anaerobes., Appl. Micro. 27, 562 (1974)
- (4) D.F. Welch, P.A. Ahlin, J.M. Matsen, Differentiation of *Haemophilus spp.* in Respiratory Isolate Cultures by an Indole Spot Test., J. Clin. Micro. 15, 216 (1982)
- (5) H.D. Isenberg, Ed., Clinical Microbiology Procedures Handbook, Vol 1., Washington, DC, ASM (1992)
- (6) B.A. Forbes, D.F. Sahm, A.S. Weissfeld, Bailey and Scott's Diagnostic Microbiology., 10th ed., St Louis, Mosby (1998)
- (7) J.F. MacFaddin, Biochemical Tests for Identification of Medical Bacteria., 3rd ed., Philadelphia, Lippincott Williams & Wilkins (2000)



Picture DMACA Indole Disks

1. *Staphylococcus aureus*
2. *Escherichia coli*

Fluka Cat. No. 405 Hippurate Disks The hippurate disks are recommended for qualitative processes of detection of organisms that have hippurate hydrolase.

The hippurate disks are recommended for qualitative detection of organisms that have hippurate hydrolase. Hippurate hydrolase promotes the hydrolysis of peptide bonds in the hippurate molecule, releasing glycine and benzoic acid as end products. The benzoic acid can be detected by using a ferric chloride indicator and it is also possible to detect glycine with ninyhydrin. However, it must be taken into account that any free amino acid will generate false positive results. The disk method is a rapid test for the presumptive identification of *Gardnerella vaginalis*, *Campylobacter jejuni*, *Listeria monocytogenes* and β -hemolytic group B *Streptococci*.

Composition

Each package contains 25, 10 mm diameter sterile filter paper disks impregnated with sodium hippurate.

Instructions

Aseptically place the hippurate disk in the brain heart infusion broth (Fluka Cat. No. 53286) inoculated with a suspect colony, incubate at 35°C for 48 hours and then separate the supernatant from the cells by centrifugation. Add 2 mL of ferric chloride reagent to 2 mL of supernatant. Shake well and check for the formation of a precipitate. If brown flocculants precipitate persists after shaking during 10 minutes, then hippurate hydrolysis can be inferred.

Preparation of ferric chloride reagent

12 g ferric chloride, 94.6 mL distilled water and 5.4 mL concentrated hydrochloric acid. Put 75 mL of distilled water in a 100 mL graduated flask. Pipette cautiously 5.4 mL of HCl into the flask and add 12 g of ferric chloride. Dissolve by warming the flask gently, swirling the contents to mix well. Bring the volume up to 100 mL with distilled water. This solution should have an orange colour.

Quality control

The table below illustrates the performance of strains used routinely for control.

Test Organisms (ATCC)	Growth	Hippurate hydrolysis
<i>Enterococcus faecalis</i> (29212)	+++	–
<i>Streptococcus agalactiae</i> (4768)	+++	+
<i>Streptococcus pyrogenes</i> (19615)	+++	–

+ = Brown flocculants precipitate persisting after shaking during 10 minutes.

– = If any visible precipitate can be dissolved by shaking.

References

- (1) S.H. Ayers, P. Rupp, Differentiation of Hemolytic *Streptococci* from Human and Bovine Sources by the Hydrolysis of Sodium Hippurate. *J. Infect. Dis.*, 30, 388 (1922)
- (2) Z. R.R. Facklam, et al., Presumptive identification of group A, B, and D *Streptococci*, *Appl. Microbiol.*, 27(1), 107 (1974)
- (3) Z. S.M. Harvy, Hippurate Hydrolysis by *Campylobacter fetus*, *J. Clin. Microbiol.* 11,435 (1980)
- (4) M. Hwang, G.M. Ederer, Rapid Hippurate Hydrolysis Method for Presumptive Identification of Group B *Streptococci*, *J. Clin. Microbiol.* 1, 114 (1975)
- (5) N.W. Luechtefeld, W.L. Wang, Hippurate hydrolysis by and triphenyltrazolium tolerance of *Campylobacter fetus*, *J. Clin. Microbiol.*, 15, 137 (1982)
- (6) P. Piot, et al., Identification of *Gardnerella (Haemophilus) vaginalis*, *J. Clin. Microbiol.*, 19 (1982)
- (7) A.E. Greenberg, R.R. Trussell, L.S. Clesceri, Eds., Standard Methods for the Examination of Water and Wastewater, 16th ed., APHA, Washington, DC (1985)
- (8) S.M. Finegold, E.J. Baron, Bailey & Scott's Diagnostic Microbiology, 8th Ed., St. Louis, MO, C.V. Mosby Co (1990)
- (9) H.D. Isenberg, Ed., Clinical Microbiology Procedures Handbook, Vol I & II, Washington, DC, ASM (1992)
- (10) P.R. Murray, et al., Manual of Clinical Microbiology, 6th ed., American Society for Microbiology, Washington D.C. (1995)
- (11) B.A. Forbes, D.F. Sahm, A.S. Weissfeld, Bailey and Scott's Diagnostic Microbiology, 10th ed., St Louis, Mosby (1998)
- (12) E.W. Koneman, S.D. Allen, W.M. Janda, P.C. Schreckenberger, W.C. Winn, Eds., Color Atlas and Textbook of Diagnostic Microbiology, 5th ed., Philadelphia: Lippincott Williams & Wilkins (1997)
- (13) Mackie and McCartney, Practical Medical Microbiology 14th ed., Vol. 2, Collee, Duguid, Fraser and Marmion, Eds., Churchill Livingstone, Edinburgh (2000)

Fluka Cat. No. 06728 Hydrogen Sulfide Test Strips (Lead Acetate Test Strips, H₂S Test Strips) Hydrogen sulfide test strips are used for detection of H₂S production by microorganisms.

A large number of bacteria can produce H₂S in small amounts from sulfur containing amino acids in carbohydrate media. When combined with lead acetate, the H₂S will produce a black precipitate, giving rise to a visible black coloured reaction on the paper strip. The lead acetate method is very sensitive, allowing the detection of trace levels of hydrogen sulphide.

Composition

Each package contains 25 sterile filter paper strips impregnated with lead acetate.

Instructions

Inoculate peptone water (Fluka Cat. No. 70179) with the suspect organism. Insert a lead acetate paper strip between the plug and inner wall of tube, above the inoculated medium and incubate at 35°C for 18-24 hours.

Quality control

Culture response after 18-24 hours at 35°C. A positive reaction appears as a blackening of the lower part of the strip. In the case of negative response, no blackening should appear.


Test Organisms (ATCC)

H₂S production

<i>Escherichia coli</i> (25922)	–
<i>Salmonella</i> serotype <i>Enteritidis</i> (13076)	+
<i>Salmonella</i> serotype <i>Typhimurium</i> (14028)	+

References

- (1) R.N. Collins, M.D. Treger, J.B. Goldsby, J.R. III Boring, D.B. Coohon, R.N. Barr. Interstate outbreak of *Salmonella newbrunswick* infection traced to powdered milk., *JAMA*, 203, 838 (1968)
- (2) J.B. Weissman, R.M.A.D. Deen, M. Williams, N. Swanston, S. Ali, An island-wide epidemic of salmonellosis in Trinidad traced to contaminated powdered milk., *West Indian Med. J.*, 26, 135 (1977)
- (3) S.M. Finegold, W.J. Martin, Bailey and Scott's Diagnostic Microbiology 6th ed., The CV. Mosby Co., St. Louis (1982)
- (4) J.F. MacFaddin, Media For Isolation-Cultivation-Identification-Maintenance of Medical Bacteria., Vol. 1, Williams and Wilkins, Baltimore (1985)
- (5) B. Rowe, N.T. Begg, D.N. Hutchinson, et al., *Salmonella* Ealing Infections Associated with Consumption of Infant Dried Milk., *Lancet*, 2, 900 (1987)



Picture Hydrogen Sulfide Test Strips

1. Control
2. *Salmonella* serotype *Typhimurium*
3. *Escherichia coli*

Fluka Cat. No. 04739 Indoxyl Strips (Acetoxyindol Strips) The indoxyls strips are a diagnostic test for the rapid detection of the acetate esterase activity of microorganisms.

The test is based on the reaction of free indoxyl groups with oxygen, which results in a colour change. Acetate esterase activity is present in species belonging to *Campylobacter* and *Branhamella catarrhalis*. The test is suitable for screening examinations and for identification of suspect colonies.

Composition

Each package contains 100 test plastic strips with an active zone saturated with 3-acetoxy indol (substrate for acetate esterase).

Storage

Store dry at +2 to +8°C. Shelf life can be extended by storing the product at -20°C.

Instructions

Wipe off several the suspect colonies grown for 18-24 hours in a Petri dish, using the paper zone of diagnostic strip. The result can be read after 3-5 minutes. For accurate results, store dry. The aluminum tube containing the strips should not be opened before the temperature is equilibrated in order to prevent condensation forming on the strips. For good test performance it is required that

there is sufficient humidity on the cultures or the cultivation media where the suspect colonies are tested. If there is insufficient humidity, the active zone of strip can be moisturised by using condensed water from a lid of the dish or by adding approximately 10 µl of distilled water. When wiping the colonies with a microbiological loop, please ensure that you do not use one of metal construction.

Interpretation of results

Negative reaction

No colour change develops at the position of wiped colony.

Positive reaction

A blue-green spot develops at the position of wiped colony.

Test Organisms (ATCC)

Escherichia coli (25922)

Branhamella catarrhalis (25238)

Result

negative

positive




Figure Indoxyl Strips

Quality control

The table above illustrates the performance of strains used routinely for control.

References

- (1) F.J. Bolton, D.R.A. Wareing, et al., Identification and Biotyping of *Campylobacters*. 29, 151. In: R.G. Board, D. Jones and F.A. Skinner. Identification Methods in Applied and Environmental Microbiology. Academic Press: London (1992)
- (2) J. Klena, A Survey of Phenotypic and Genetic Methods Used to Identify and Differentiate Thermotolerant *Campylobacter spp.* Strains, report of Ministry of Health Department of Plant and Microbial Sciences University of Canterbury (2001)

Fluka Cat. No. 78719 Kovac's Reagent Strips (Indole Reagent Test strips according to Kovac) The Kovac's reagent strips are used for determination of the ability of microorganisms, primarily *Enterobacteriaceae*, to use tryptophanase.

Tryptophanase, present in species such as *Escherichia coli*, cleaves tryptophan to indole and α -aminopropionic acid. The p-aminobenzaldehyde present in the reagent reacts with the indole to form a pink complex.

Composition

Each package contains 25 sterile filter paper strips impregnated with Kovac's reagent. The Kovac's reagent is prepared by dissolving 10 g of p-aminobenzaldehyde in 150 mL of isoamylalcohol and then slowly adding 50 mL of concentrated hydrochloric acid.

Instructions

Indole production by the organism is observed by inserting the Kovac's reagent strips between the plug and inner wall of the tube, above the inoculated peptone water (Fluka Cat. No. 70179) and incubating at 35°C for 18-24 hours.

Quality control

Culture response appears after incubation in peptone water (Fluka Cat. No. 70179) for 18-24 hours at 35°C. Pink coloration on the lower part of the strip should be considered a positive reaction, as the microorganisms have tryptophanase activity. For negative reactions there is no colour change.

References

- (1) J.F. MacFaddin, *Biochemical Tests for Identification of Medical Bacteria*, 2nd ed., Williams and Wilkins, Baltimore (1980)
- (2) A.E. Greenberg, R.R. Trussell, L.S. Clesceri (Eds.), *Standard Methods for the Examination of Water and Wastewater*, 16th ed., A.P.H.A, Washington, D.C. (1985)

Test Organisms (ATCC)


Escherichia coli (25922)

Enterobacter aerogenes (13048)

Indole production

+

-



Picture Kovac's Reagent Strips

1. Control
2. *Escherichia coli*
3. *Staphylococcus aureus*

Fluka Cat. No. 56348 β-Lactamase Strips The β-Lactamase strips are suitable for the rapid acidimetric detection of β-lactamase activity in microorganisms.

The β-lactamase test is based on hydrolysis of the β-lactam ring in benzylpenicillin, which results in the production of penicilloic acid. This process causes acidification of the bacterial suspension, therefore changing the colour of the acid-base indicator. The acid-base test for β-lactamase activity is suitable only for the detection of *Haemophilus influenzae*, *Neisseria gonorrhoeae* and *Staphylococcus spp.* This test is not suitable for the detection of the β-lactamase activity in microorganisms such as *Branhamella catarrhalis*, *Enterococcus faecalis*, *Neisseria meningitis*, *Enterococcus spp.* and others.

Composition

Each package contains 100 test strips with an active zone saturated with benzylpenicillin and an acid-base indicator.

Storage

Store dry at +2 to +8°C. The shelf life of this product can be extended by storing it at –20°C.

Instructions

Solid Media Test

Wipe off several suspect colonies from a Petri dish by the function zone of a diagnostic strip. Mark the strip and incubate at room temperature. The results can be read after 2-10 minutes. Test performance requires that there is sufficient humidity on the cultures or the

cultivation media where the suspect colonies are tested. If sufficient moisture is not present, then the active zone of the strip can be moisturised by using condensed water condensed on the lid of the dish or by adding approximately 10 µl of distilled water.

Liquid Media Test

Prepare approximately 0.5-1 mL of bacterial suspension in a saline solution (2-4 loops). Insert the test strip in the test tube with the prepared bacterial suspension, shake and incubate at room temperature. Read result after 2-10 minutes.

Interpretation of results

Negative reaction

Solution remains red or no colour change develops at the position of wiped colony.

Positive reaction

Solution turns yellow or a blue-green spot develops at the position of wiped colony.

References

- (1) R. Bonnet, C. Chanal, E. Ageron et al., Inducible AmpC β-Lactamase of a New Member *Enterobacteriaceae*, *Antimicrob. Agents Chemother.*, 46, 3316 (2002)
- (2) L. Shan, et al., Kinetic Analysis of an Inhibitor-Resistant Variant of the OHIO-1 β-lactamase, an SHV-Family Class A Enzyme, *Biochem. J.*, 333, 395, Printed in Great Britain (1998)




Figure β-Lactamase Strips

Fluka Cat. No. 08086 Nitrate Reagent Disks The nitrate reagent disks are used to detect the reduction of nitrate.

The nitrate test involves the detection of the enzyme nitrate reductase, that causes the reduction of nitrate to nitrite in the presence of a suitable electron donor. Nitrite can be tested by using an appropriate colorimetric reagent.

Composition

Each package contains 50, 6mm diameter sterile filter paper disks impregnated with nitrate reagent.

Storage

Store in the freezer below 4°C in the containers provided with the product. Allow to equilibrate to room temperature before opening. Return to freezer storage immediately after use.

Instructions

The test culture should be grown on suitable agar medium plate containing the nitrate substrate. Place the nitrate reagent disks on the suspected colony.

Principle

Reduction of nitrate (NO_3) to nitrite (NO_2) and to nitrogen gas (N_2) usually takes place under anaerobic conditions (1). Most facultative anaerobes can reduce nitrate in the absence of oxygen and almost all *Enterobacteriaceae* are able to reduce nitrate. This anaerobic respiration is an oxidation process in which inorganic substances provide oxygen to serve as an electron acceptor and provide energy (2). Depending on the bacterial species, the nitrate reduction results in the production of various end products, the most common of which is molecular nitrogen by way of nitrite reduction (2). Depending upon the environmental conditions, these products are usually not further oxidised or assimilated into cellular metabolism, but are excreted into the surrounding medium.

Interpretation of results

Colour change to red-pink of the disk indicates positive nitrate reduction reaction.

Quality control

The table below illustrates the performance of strains used routinely for control.

Test Organisms (ATCC)

Test Organisms (ATCC)	Growth	Nitrate reduction
<i>Acinetobacter calcoaceticus</i> (19606)	+++	-
<i>Enterobacter aerogenes</i> (13048)	+++	+
<i>Escherichia coli</i> (25922)	+++	+
<i>Salmonella</i> serotype Typhimurium (14028)	+++	+

References

- (1) Jr. M.J. Pelczar, R.D. Reid, Microbiology, 2nd ed., McGraw-Hill, New York, 567 (1965)
- (2) R.Y. Stainer, M. Douderoff, E.A. Adelberg, The Microbial World, 2nd ed., Prentice-Hall, 116-117 (1963)
- (3) Wideman P.A., Citronbaum D.M., Sutter V.L., Simple Disk Technique for Detection of Nitrate Reduction in Anaerobic Bacteria, J Clin Micro, 5(3):315-9 (1977)
- (4) Lennette E.H., Balows A., Hausler W.J., et al., Manual of Clinical Microbiology, 4th ed. Washington, DC: ASM (1985)
- (5) Baron E.J., Finegold S.M., Bailey and Scott's Diagnostic Microbiology, 8th ed., St Louis: Mosby, (1990)



Picture Nitrate Reagent Disks

1. *Acinetobacter calcoaceticus*
2. *Salmonella typhimurium*

Fluka Cat. No. 49862 Nitrocefin Disks For the rapid detection of β -lactamase enzymes in isolated colonies of *Neisseria gonorrhoeae*, *Moraxella catarrhalis*, *Staphylococcus spp.*, *Haemophilus influenzae* and anaerobic bacteria.

Composition

Each package contains 50, 6mm diameter filter paper disks impregnated with nitrocefin in a light resistant plastic vial.

Storage

Store in the freezer below -10°C in the containers provided. Allow to equilibrate to room temperature before opening and return to freezer storage immediately after use.

Instructions

Place the required number of nitrocefin disks into a clean, empty Petri dish or onto a microscope slide. Disks may be moistened with one drop of deionised water, but take care not to over-moisten. Using a sterile loop or an applicator stick, remove several similar well-isolated colonies and smear them onto the surface of a disk. Alternatively, moisten the disk with one drop of deionised water and then, holding the disk with forceps, wipe across a colony on an agar plate. Observe the inoculated disk for the development of a red colour.

Interpretation of results

Positive

Development of a red colour in the area of the disk where the culture was applied. Note the colour change does not normally develop over the whole of the disk.

Negative

No colour change.

A positive result should be interpreted as resistance to penicillin or cephalosporin activity. The susceptibility should be confirmed by standard growth-dependent susceptibility testing methods. A negative result implies but does not guarantee susceptibility.

Quality control

The table below illustrates the performance of strains used routinely for control.

User quality control

Check for signs of deterioration. Quality control must be performed with at least one organism to demonstrate a positive reaction and at least one organism to demonstrate a negative reaction. Do not use the product if the reactions with the control organisms are incorrect.

Test Organisms (ATCC)

Test Organisms (ATCC)	Result
<i>Haemophilus influenzae</i> (35036)	negative
<i>Neisseria gonorrhoeae</i> (31426)	positive
<i>Staphylococcus aureus</i> (11632)	positive
<i>Escherichia coli</i> (25922)	negative

Limitations

It is recommended that biochemical and/or serological tests are performed on colonies from pure cultures to confirm identification.

For most bacterial strains a positive results will develop within 5 minutes. However, a positive reaction for some *Staphylococci* and anaerobic species may take up to 60 minutes to develop.

Detection of staphylococcal β -lactamase is enhanced by testing growth from a medium containing sub-inhibitory concentrations of a β -lactam antibiotic.

References

- (1) MacFaddin J.F., *Biochemical Tests for the Identification of Medical Bacteria.*, 3rd ed. Philadelphia: Lippincott Williams & Wilkins (2000)
- (2) Murray P.R., Baron E., Pfaller M., Tenover F., Tenover R., *Manual of Clinical Microbiology*, 7th ed. Washington, DC: ASM, (1999)
- (3) Tu K.K., Jorgensen J.H., Stratton C.W., A Rapid Paper-Disk Test for Penicillinase., *Am J. Clin. Pathol.*, 75:557-9 (1981)
- (4) Escamilla J., Susceptibility of *Haemophilus influenzae* to Ampicillin as Determined by Use of a Modified One-Minute β -Lactamase Test., *Antimicrob. Ag. Chemo.*, 9:196-8 (1976)
- (5) Thornsberry C., Biddle J.W., Kirven L.A., Penicillin Resistance in *Neisseria gonorrhoeae* due to β -Lactamase Production., *Microbios* 20:39-46 (1977)

Fluka Cat. No. 49940 ONPG Disks (2-Nitrophenyl β -D-galactopyranoside Disks, β -Galactosidase Test Disks) ONPG disks are used to detect the presence of β -galactosidase, an enzyme found in lactose-fermenting organisms.

Lactose utilization depends upon two enzymes: a β -galactoside permease, not present in late lactose fermenters, which catalyzes transport of lactose into the cell and β -galactosidase, which breaks down lactose into galactose and glucose. β -Galactosidase is not lactose specific and can act on simple galactosides including the ONPG (o-nitrophenyl- β -D-galactopyranose) substrate. ONPG hydrolysis results in the release of galactose, and o-nitrophenol, a yellow chromogen. The test substrate ONPG does not depend on an induced or constitutive permease enzyme to enter the cell. Therefore, the reactions are rapid and occur within a 24-hour period even for late lactose fermenters. To group *Enterobacteriaceae* the ability of fermenting lactose is routinely used.

Composition

Each package contains 50, 6 mm diameter sterile filter paper disks impregnated with o-nitrophenyl- β -D-galactopyranose.

Instructions

Place one ONPG disk into a sterile test tube. Add 0.1 mL of sterile 0.85% (w/v) sodium chloride solution (physiological saline). Pick up the colony to be tested with a sterile loop and emulsify it in the tube containing the disk and physiological saline. Incubate at 35°C. To detect active lactose fermenters observe the tube at hourly intervals up to 6 hours. To detect late lactose fermenters, incubate the negative tubes for up to 24 hours.

Quality control


Culture characteristics after growth in a 0.85% (w/v) sodium chloride solution with an ONPG disk during 4 hours at 35°C.

References

- (1) W.L. Gaby, C. Hadley, J. Bact., 74, 356 (1957) J. Sanbrook, E.F. Fritsch, T. Maniatis, Molecular Cloning: A Laboratory Manual 2nd ed., Cold Spring Harbor, NY (1989)
- (2) S.R. Maloy, J.E. Conran (Jr.), D. Freifelder, Microbial Genetics 2nd ed. Jones and Bartlett Boston, MA (1994)
- (3) V.E. Becker, H.J. Evans, The influence of Monovalent Cations and Hydrostatic Pressure on β -Galactosidase Activity, Biochim. Biophys. Acta, 191, 95 (1969)
- (4) M.C. Neville, G.N. Ling, Synergistic Activation of β -Galactosidase by Na⁺ and Cs⁺, Arch. Biochem. Biophys., 118, 596, (1967)
- (5) J. Lederberg, The β -Galactosidase of *Escherichia coli*, strain K-12., J. Bact., 60, 381 (1950)

Picture ONPG Disks (from left to right)

1. Positive Colony
2. Negative Colony
3. Negative Colony
4. Negative Control



Test Organisms (ATCC)

ONPG Hydrolysis

<i>Citrobacter freundii</i> (8090)	+
<i>Enterobacter aerogenes</i> (13048)	+
<i>Escherichia coli</i> (25922)	+
<i>Proteus vulgaris</i> (8427)	-
<i>Salmonella arizonae</i> (13314)	+
<i>Salmonella serotype Typhimurium</i> (14028)	-

Fluka Cat. No. 74042 Optochin Disks The optochin test is a useful diagnostic tool for identification/differentiation of *Pneumococci* and *viridans Streptococci*.

Optochin (ethyl hydrocuprein hydrochloride) is inhibitory for pneumococcal growth whereas other *Streptococci* show good growth or a very small zone of inhibition. Further tests are required for the diagnosis of pneumococcal infections, since α -haemolytic (*viridans*) *Streptococci* and *Pneumococci* (*Streptococcus pneumoniae*) cannot be easily differentiated on a blood agar plate (show both a partial clearing of blood and a greenish discolouration; α -haemolysis). The optochin test disks are suitable for the examination of specimens like sputum, pleural fluid, lung aspirate, urine or blood. The correlation between bile solubility and full optochin susceptibility for the differentiation of *Streptococcus pneumoniae* from other *Streptococci* was shown by Bowers and Jeffries (1).

Composition

Each package contains 50 sterile filter paper disks impregnated with optochin.

Instructions

Gram staining test (Fluka Cat. No. 77730) is required, followed by the optochin sensitivity test. Prepare tryptone soya agar (Fluka Cat. No. 22091) with blood or blood agar (Fluka Cat. No. 70133) plates and streak pure culture of organism to be tested across one half of the plate. Streak a known *Pneumococcus* culture across the other half of the plate as positive control and immediately place the optochin disks in the centre of the two halves of the plate and incubate at 37°C. Observe for zone of inhibition around the disks.

Quality control

Culture response after 18-24 hours at 35-37°C on seeded tryptone soya agar (Fluka Cat. No. 22091) with blood using optochin disks.

Test Organisms (ATCC)

Diameter of zone of inhibition

<i>Streptococcus pneumoniae</i> (6303)	15 mm
<i>Streptococcus pyrogenes</i> (19615)	13 mm

References

- (1) E.F. Bowers, L.R. Jeffries, J. Clin. Path., 8, 58 (1995)
- (2) Bouvet A., Grimont F., Grimont P.A.D., *Streptococcus defectivus* sp. nov. and *Streptococcus adjacens* sp. nov., Nutritionally Variant *Streptococci* from Human Clinical Specimens., Int. J. Syst. Bacteriol., 39, 290 (1989)
- (3) Baron E.J., Peterson L.R., Finegold S.M., Bailey and Scott's Diagnostic Microbiology., 9th edition. St. Louis, Mosby (1994)

Fluka Cat. No. 40560 Oxidase Strips The oxidase strips test is a diagnostic test for the detection of the cytochrome oxidase activity in microorganisms within 1 minute.

The oxidase strips test is an important differential procedure which should be performed on all Gram-negative bacteria that are to be identified. The cytochrome oxidase present in most Gram-negative bacteria triggers the reaction of N,N-dimethyl-p-phenylenediamine with α -naphthol, forming indophenol blue.

Composition

The kit contains 100 plastic strips with a paper zone saturated with a solution of N,N-dimethyl-1,4-phenylenediamine and α -naphthol.

Storage

Store dry at +2 to +8°C. The shelf life of this product can be extended by storing it at -20°C.

Instructions

Wipe off several suspect colonies from a Petri dish by the paper zone of diagnostic strip. Read result after 1 minute.

For accurate results, store dry. The aluminum tube containing the strips cannot be opened before the temperature has equilibrated to room temperature to prevent condensation of air humidity onto the strips. However, test performance requires that there is sufficient humidity on the cultures or the cultivation media where the suspect colonies are tested. If sufficient moisture is not available, the active zone of strip can be moisturised by using water condensed on to the lid of the dish or by adding approximately 10 µl of distilled water. When wiping the colonies with a microbiological loop, please ensure that you do not use one of metal construction.

Principle and Interpretation of Results

Gordon and McLeod (1) introduced the oxidase test for identifying *Gonococci* based upon the ability of certain bacteria to produce indophenol blue from the oxidation of dimethyl-p-phenylenediamine and α -naphthol. Later, Gaby and Hadley (2) introduced a more sensitive method by using N,N-dimethyl-p-phenylenediamine oxalate, establishing that *Staphylococci* are Gram-negative. The cytochrome oxidase present in most Gram-negative bacteria triggers the reaction of N,N-dimethyl-p-phenylenediamine with α -naphthol, forming indophenol blue.

Oxidase test is mainly used to differentiate

1. Oxidase positive *Neisseria* from other Gram-negative diplococci
2. Oxidase positive *Aeromonas hydrophila* from *Escherichia coli* (Gram-negative)
3. Oxidase positive *Plesiomonas shigelloids* from *Shigella sonnei* (Gram-negative)

Cytochrome oxidase production may be inhibited by acid production and false negative reaction may be given by *Vibrio*, *Aeromonas*, and *Plesimonas spp* when grown on a medium containing fermentable carbohydrate such as MacConkey Agar (Fluka Cat. No. 70143). Colonies taken from media containing nitrate may give unreliable results. The loss of activity of the oxidase reagent is caused by auto oxidation which may be avoided by adding 0.1% ascorbic acid (Fluka Cat. No. 95209).

Negative reaction

No colour change develops at the position of wiped colony.

Positive reaction

A dark blue or black spot develops at the position of wiped colony.

Quality control

The list below illustrates control strains in routine use.

Test Organisms (ATCC)

Escherichia coli (25922)

Pseudomonas aeruginosa (27853)

Result

negative

positive




Figure Oxidase Strips

References

- (1) J.Gordon, J.W. McLeod, J. Path.Bact., 31, 185 (1928)
- (2) W.L. Gaby, C. Hadley, J. Bact., 74, 365 (1957)
- (3) K.J. Steel, J. Appl. Bact., 25, 445 (1962)

Fluka Cat. No. 70439 Oxidase Test Disks The oxidase test is an important differential procedure which should be performed on all Gram-negative bacteria that are to be identified.

Composition

Each package contains 50 disks impregnated with N,N-dimethyl-p-phenylenediamine oxalate and α -naphthol.

Storage

Store dry at +2 to +8°C. The shelf life of this product can be extended by storing it at -20°C.

Instructions

With an inoculating loop or a toothpick, touch and spread a well isolated colony on an oxidase disk. The reaction is observed within 2 minutes at 25-30°C.

Principle and Interpretation of Results

Gordon and McLeod (1) introduced the oxidase test for identifying *Gonococci* based upon the ability of certain bacteria to produce indophenol blue from the oxidation of dimethyl-p-phenylenediamine and α -naphthol. Later, Gaby and Hadley (2)

introduced a more sensitive method by using N,N-dimethyl-p-phenylenediamine oxalate, establishing that *Staphylococci* are Gram-negative. The cytochrome oxidase present in most Gram-negative bacteria triggers the reaction of N,N-dimethyl-p-phenylenediamine with α -naphthol, forming indophenol blue.

Oxidase test is mainly used to differentiate

1. Oxidase positive *Neisseria* from other Gram-negative diplococci.
2. Oxidase positive *Aeromonas hydrophila* from *Escherichia coli* (Gram-negative)
3. Oxidase positive *Plesiomonas shigelloids* from *Shigella sonnei* (Gram-negative)

Cytochrome oxidase production may be inhibited by acid production and false negative reaction may be given by *Vibrio*, *Aeromonas*, and *Plesimonas* species when grown on a medium containing fermentable carbohydrates such as MacConkey Agar (Fluka Cat. No. 70143). Colonies taken from media containing nitrate may give unreliable results. The loss of activity of the oxidase reagent is caused by auto oxidation which may be avoided by adding 0.1% ascorbic acid (Fluka Cat. No. 95209).


Reaction takes place within 2 minutes at 25-35°C.

Organisms (ATCC)	Reaction	Colour
<i>Pseudomonas aeruginosa</i> (27853)	positive	deep purple blue
<i>Staphylococcus aureus</i> (25923)	negative	-
<i>Neisseria gonorrhoeae</i> (19424)	positive	deep purple blue
<i>Escherichia coli</i> (25922)	negative	-

References

- (1) J.Gordon, J.W. McLeod, J. Path.Bact., 31, 185 (1928)
- (2) W.L. Gaby, C. Hadley, J. Bact., 74, 365 (1957)
- (3) K.J. Steel, J. Appl. Bact., 25, 445 (1962)

Picture Oxidase Disks



Fluka Cat. No. 67886 PYRase Strips (Pyrrolidonyl Peptidase Strips)

Diagnostic test for the rapid differentiation of *Enterococci* from the group D *Streptococci* and differentiation of *Streptococcus pyogenes* from other haemolytic *Streptococci*.

Composition

The kit contains 50 plastic strips with a paper zone saturated with chromogenic substrate for the detection of pyrrolidonyl peptidase and 1.0 mL of developing reagent.

Storage

Store dry at +2 to +8°C. The shelf life of this product can be extended by storing it at -20°C.

Instructions

Wipe off several suspect colonies from a Petri dish by the paper zone of diagnostic strip. After one minute, apply approximately 10 µl of developing reagent to the paper zone using a plastic loop or micropipette. The reaction is observed within 1 minute.

Interpretation of results

Negative reaction

No colour change develops at the position of wiped colony.

Positive reaction

A red spot develops at the position of wiped colony.

Use PYRase test due to the recommendation of NRL for *Streptococci* and *Enterococci*.

Positive (change to red colour)

Enterococcus
Streptococcus pyogenes
Lactococcus lactis
Aerococcus
Gamella (most strains)

Negative (no colour change)

Streptococcus bovis
Streptococcus equinus
Beta Minute haemolytic *Streptococci*:
Streptococcus anginosus
Streptococcus intermedius
Streptococcus constellatus
Other *Lactococcus* species
Leuconostoc
Pediococcus
Lactobacillus (sporadic positive test)
Other species of *viridans Streptococci*




Figure PYRase Strips

References

- (1) Godsey et al., Abstr. ASM Annu. Meet., C84, 276 (1981)
- (2) R.R. Facklam, L.G. Thacker, B. Fox, L. Eriquez, Presumptive Identification of *Streptococci* with a New Test System, J. Clin. Micro., 15, 987 (1982)
- (3) B.L. Wasilauskas, K.D. Hampton, Evaluation of the Strep-A-Fluor Identification Method for Strep A *Streptococci*, J. Clin. Micro. 20, 1205 (1984)
- (4) D.M. Yajko, J. Lawrence, P. Nassos, J. Young, K.W. Hadley, Clinical trial comparing bacitracin with Strep-A-Chek for accuracy and turnaround time in the presumptive identification of *Streptococcus pyogenes*, J. Clin. Micro., 24, 431 (1986)
- (5) L.P. Gordon, M.A.S. Damm, J.D. Anderson, Rapid Presumptive Identification of *Streptococci* Directly from Blood Cultures by Serologic Tests and the L-Pyrrolidonyl-β-Naphthylamide Reaction, J. Clin. Micro. 25, 238 (1987)
- (6) P.R. Murray, E.J. Baron, M.A. Tenover, F.C. Tenover, R.H. Tenover, Eds. Manual of Clinical Microbiology, 7th ed., Washington, D.C., ASM Press (1999)
- (7) J.F. MacFaddin, Biochemical tests for the Identification of Medical Bacteria, 3rd ed., Philadelphia, Lippincott Williams & Wilkins (2000)

Fluka Cat. No. 74041 Sterility Indicator (Steam Sterilization) These indicators consist of ca 1 million *Bacillus stearothermophilus* (ATCC 7953) spores impregnated on paper strips, individually placed into envelopes.

Bacillus stearothermophilus is a thermophilic species that can grow at temperatures such as 65°C and higher. Its spores are an excellent tool to monitor autoclave performance, as they are highly resistant to temperature. These indicators are similar to those specified by the United States military specification MIL-S-36586 are GMP requirements of the United States FDA.

Composition

Each package contains 25 test strips. One paper strip is impregnated with 1 million spores of *Bacillus stearothermophilus* (ATCC 7953).

Instructions

Place the indicators in such a manner that they are exposed to the steam in representative places of the media and the equipment. For routine evaluation of each radiated lot, a standard procedure should be established. After steam sterilization open the envelope with the strip inside by using strict aseptic techniques. Inoculate the strip after autoclaving in a tube with sterile CASO broth (Fluka Cat. No. 22098) and incubate at 55-60°C up to 7 days. Also, an unexposed spore strip should be inoculated at the same time in another tube with sterile CASO broth.

Quality control

Culture response observed after 7 days at 55-60°C in CASO Broth (Fluka Cat. No. 22098):

Spore Strip	Result
Exposed to steam	no growth
Unexposed Luxuriant	growth

Fluka Cat. No 74146 Sterile Disks The sterile disks can be used to test a variety of antibiotics, carbohydrates, substrates, antiseptics on bacteria in Petri dishes.



The sterile disks can be used to test a variety of antibiotics, carbohydrates, substrates and antiseptics on bacteria in Petri dishes. Soak a disk in a solution or apply some solution on the disks. Allow to dry and place it in an inoculated agar plate. Each disk will absorb exactly the same amount of liquid.

Picture Sterile Disks

Fluka Cat. No. 05290 Sterility Indicator (Radiation Sterilization)

These indicators consists of ca 1 million *Bacillus pumilus* (ATCC 27142) spores impregnated on paper strips, individually placed into envelopes.

Bacillus pumilus was chosen due to its resistance to radiation. Its spores are an excellent tool to monitor the efficiency of radiation sterilization since they are highly resistant. The above mentioned indicators are similar to those specified by the U.S. military specification MIL-S-36586 and are GMP requirements of the United States FDA.

Composition

Each package contains 25 test strips. One paper strip is impregnated with 1 million spores from *Bacillus pumilus* (ATCC 27142).

Instructions

Place indicators, in such a manner that they are exposed to the same radiation as the media and the equipment. For the routine evaluation of each radiated lot, a standard procedure should be established. After radiation sterilization, open the envelope with the strip using rigid aseptic techniques. Inoculate strip in a tube with sterile CASO broth (Fluka Cat. No. 22098) and incubate at 35-37°C up to 7 days. Also, an unexposed spore strip should be inoculated at the same time in another tube with sterile CASO broth.

Quality control

Culture response observed after 7 days at 35°C in CASO Broth (Fluka Cat. No. 22098):

Spore Strip	Result
Exposed to radiation	No growth
Unexposed	Luxuriant growth



Picture Sterility Indicator

Fluka Cat. No. 75744 Tributyrin-Strips (TRIBU Strips) Tributyrin strips are a diagnostic test for the differentiation of *Branhamella* and *Neisseria*. The test principle is the enzymatic hydrolysis of tributyrin, a reaction that causes a colour change of acid-base indicator. The results can be read after 18-20 hours.

Composition

Each package contains 300 test strips saturated with tributyrin and acid-base indicator.

Storage

Store dry at +2 to +8°C. The shelf life of this product can be extended by storing it at –20°C.

The strips are delivered sterilized. Observe the maintenance of sterility when the strips are used repeatedly. To obtain accurate results, avoid moisture during storage. Allow to equilibrate to room temperature before opening the container. When this is not observed, the product can become moistened by condensation, and will deteriorate.

Instructions

Using a sterile forceps, throw one tributyrin strip into a 1mL suspension of the test strain in buffered saline (pH 7.2). Incubate the test sample at 37°C (without CO₂). Preliminary results can be read after several hours when the red colour changes to yellow (positive result). The final result can be obtained after 18-20 hours incubation.

Interpretation of results

Negative reaction

Red colour did not change to yellow (*Neisseria*)

Positive reaction

Red colour change to yellow (*Branhamella*)

Identification diagram for *Branhamella*

Tributyrim Reduction of nitrates	
+	–
<i>Branhamella</i>	<i>Neisseria</i>
saccharides +/- β-lactamase +/-	saccharides +/-
+	–
adapted parasite	opportunistic pathogen

Quality control

The list below illustrates control strains in routine use.

Test Organisms (ATCC)	Result
<i>Neisseria gonorrhoeae</i> (19424)	negative
<i>Branhamella catarrhalis</i> (25238)	positive

References

- Berger U., Über die Spaltung von tributyrin durch *Neisseria*., Arch Hyg. Bakteriol., 146:388-391 (1962)
- Kuzmanská P., Laboratorní průkaz kokobacilli a Gramnegativních koků, Avicenum Praha, 52-57, (1987)
- Janda W.M. and P. Ruther. B. CAT CONFIRM; A rapid test for confirmation of *Branhamella catarrhalis*. J. Clin. Microbiol., 27:1390-1391 (1989)
- August M.J., et al., Cumitech 3A; Quality Control and Quality Assurance Practices in Clinical Microbiology, Coordinating ed., A.S. Weissfeld. American Society for Microbiology, Washington D.C. (1990).
- Perez J.L., et al. Butyrate Esterase (Tributyrim) Spot Test, a Simple Method for Immediate Identification of *Moraxella (Branhamella) catarrhalis*. J. Clin. Microbiol., 28: 2347-2348 (1990)
- Murray P.R., et al. Manual of Clinical Microbiology, 6th ed. American Society for Microbiology, Washington D.C. (1995)
- Koneman E.W., et al. Color Atlas and Textbook of Diagnostic Microbiology, 5th ed. J.B. Lippincott Company, Philadelphia, PA, (1997)
- Forbes B.A., et al. Bailey and Scott's Diagnostic Microbiology, 10th ed. C.V. Mosby Company, St. Louis, MO, (1998)

Fluka Cat. No. 77148 Differentiation Disks X Factor / Fluka Cat. No. 89788 V Factor / Fluka Cat. No. 08482 X + V Factors Use for the presumptive identification of *Haemophilus* species on the basis of their requirements for X or V growth factor or both.

The members of the genus *Haemophilus* require hemin (X factor) and/or nicotinamid-adenin-dinucleotid (V-factor) to growth. The need for either one or both factors provides a way for differentiation of these organisms.

Composition

Each package contains 50 test disks. Sterile 6mm diameter filter paper disks impregnated with:

- hemin and nicotinamide adenine dinucleotide (X + V factor disks; Fluka Cat. No. 08482).
- hemin (X factor disks; Fluka Cat. No. 77148).
- nicotinamide adenine dinucleotide (V factor disks; Fluka Cat. No. 89788).

Instructions

Inoculate the surface of a blood agar (Fluka Cat. No. 70133) or brain heart infusion agar (Fluka Cat. No. 70138) plate with the test organisms either by streaking or surface spreading. Aseptically, place the X and V factor disks on the plate. Incubate the plates at 35-37°C for 24-48 hours.

Recommended Disk Positions on the Agar Plate

Disk	Place
X factor disk	12 o'clock
V factor disk	4 o'clock
X + V factor disk	8 o'clock

Observe the growth in the neighbourhood of the disk. The test organism requiring X factor grows only in the vicinities of X disks. Those who require V factor grows only in the vicinities of V disks. The ones who require V and X factor grows only in the vicinities of X + V disks.

Note: Use known strains of *Haemophilus influenza* to monitor the performance of the differentiation disks and the medium. Do not use an heavy suspension of the test organisms as X- or V-factor carry over from the primary growth medium may take place.

Quality control

Culture response observed on brain blood agar (Fluka Cat. No. 70133) plate or brain heart infusion agar (Fluka Cat. No. 70138) after 24-48 hours at 35-37°C.

Test Organisms (ATCC)	Without growth factor	X factor	Y factor	X +Y factor
<i>Haemophilus influenzae</i> (35056)	–	–	–	+
<i>Haemophilus parainfluenza</i> (7901)	–	–	+	+
<i>Haemophilus ducreyi</i> (27722)	–	+	–	+
<i>Bordetella pertussis</i> (13048)	+	+	+	+

References

- (1) D.J. Davis, J. Infect. Dis., 21, 392 (1917)
- (2) T. Thjotta, O.T. Avery, J. Exp. Med., 34, 97 (1921)
- (3) A. Lwoff, M. Lwoff, Proc. R. Soc., 122, 352 London (1937)
- (4) R.H. Parker, P.D. Hoepflich, Am. J. Clin. Pathol., 37, 319 (1962)
- (5) E.L. Biberstein, P.D. Mini, M.G. Bills, J. Bacteriol., 86, 814 (1963)
- (6) D.C. White, S. Granick, J. Bacteriol., 85, 842 (1963)
- (7) P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, R.H. Tenover, Manual of Clinical Microbiology, 6th ed. ASM, Washington, D.C. (1995)
- (8) A.E. Greenberg, R.R. Trussell, L.S. Clesceri (Eds.), Standard Methods for the Examination of Water and Wastewater, 16th ed., A.P.H.A, Washington, D.C. (1985)
- (9) M. Kilian, *Haemophilus*. In Manual of Clinical Microbiology, Edited by A. Balows, W.J. Hausler, K.L. Herrmann, H.D. Isenberg, H.J. Shadomy., p 463. Washington, DC: American Society for Microbiology (1991)

Antimicrobial Susceptibility Disks recommended by the WHO Expert Committee on Biological Standardization.

Antimicrobial susceptibility disks are paper disks impregnated with antimicrobial substances, used for identification of suspect organisms or for the study of resistance to antibiotics. The *in vitro* susceptibility of pathogenic organisms is measured by using techniques such as disk-agar diffusion or disk broth elution. For disk agar diffusion, the bacterial susceptibility is ascertained by measuring the zone of bacterial inhibition around the disks on a agar surface of media such as Mueller Hinton agar (Fluka Cat. No. 70191). Disk-broth elution is associated with an automated rapid susceptibility test system and employs fluid media, such as the Mueller Hinton broth (Fluka Cat. No. 70192). When the disk is placed into the medium, the antimicrobial substance is eluted. The resulting changes in the bacterial growth are then measured by using a photometer.

The disks are manufactured under aseptic conditions by a dry process, resulting in an extra stable material (3 to 4 years, at room temperature). For better identification, each microbial agent disk is marked with a different colour.

The quality control is made using stock cultures of *Staphylococcus aureus* (ATCC 25923), *Escherichia coli* (ATCC 25922), *Pseudomonas aeruginosa* (ATCC 27853)

and *Staphylococcus faecalis* (ATCC 29212), obtained from reliable sources. Factors that seriously affect the accuracy and precision of the disk susceptibility assay must be rigorously controlled. These include: preparation of the media, thickness of the agar, inoculum density, or any change in test procedure, such as using a new lot of disks, different medium or even staff.

The susceptibility disks can be also used for testing the presence of broad-spectrum β -lactamase, for example in species belonging to the *Enterobacteriaceae*. β -lactamases hydrolyze penicillins and oxyminocephalosporins (cefotaxime, ceftazidime and ceftriaxone), but are susceptible to certain inhibitors such as clavulanic acid [1-2]. Their presence can be detected by performing a synergy test, that combines sensitivity with susceptibility in one single plate. The sensitivity test is done by using Mueller-Hinton medium. For the susceptibility test, ceftazidime, cefotaxime and amoxicillin-clavulanic acid disks are placed at a distance of 1.5 cm from each other. Incubate overnight at 37 °C. The synergic action of all of these factor makes the inhibition areas not to be perfect circles, as in a normal susceptibility test. In this case, the enlargement of the ceftazidime and cefotaxime inhibition areas towards the disk containing amoxycillin and clavulanic acid is an indication of the presence of broad spectrum β -lactamases [3].

References

- (1) Brown R.P., Aplin RT, Schofield C.J., Inhibition of TEM-2 β -lactamase from *Escherichia coli* by Clavulanic Acid: Observation of Intermediates by Electrospray Ionization Mass Spectrometry., *Biochemistry*. 1996 Sep 24;35(38):12421-32.
- (2) Philippon L.N., Naas T., Bouthors A.T., Barakett V., Nordmann P., OXA-18, a Class D Clavulanic Acid-Inhibited Extended-Spectrum β -lactamase from *Pseudomonas aeruginosa*., *Antimicrob Agents Chemother*. 1997 Oct;41(10):2188-95.
- (3) R. Skov, N. Frimodt-Møller, F. Espersen, *In vitro* Susceptibility of *Staphylococcus aureus*, *APMIS* 2002 110:7-8 p. 559



Picture Susceptibility test to of *Staphylococcus ssp.* isolated from food grown on Mueller Hinton agar (Fluka Cat. No 70191). Disk color code: ampicillin: red/blue; cefalexin: blue/beige; norfloxacin: pink/gray; erythromycin: green; ciprofloxacin: yellow/gray; imipenem: orange/blue.

Fluka Cat. No.	Antibiotic	Disk content	Colour	Acc. NCCLS
42456	Amikacin	30 µg	yellow/green	✓
16527	Amoxicillin	25 µg	blue/white	
1601	Amoxicillin + Clavulonic acid	20 + 10 µg	blue/white/orange	✓
8541	Ampicillin	10 µg	red/blue	✓
68601	Azithromycin	15 µg	green/orange	✓
52856	Cefachlor	30µg	gray/pink	✓
68611	Cefalexin	30 µg	blue/beige	
30321	Cefotaxim	30 µg	blue/dark green	
89867	Ceftazidime	30 µg	beige/yellow	✓
90424	Ceftriaxone	30 µg	pink/dark blue	✓
92241	Cefuroxime	30 µg	white/blue	✓
7651	Chloramphenicol	30 µg	white	✓
8587	Ciprofloxacin	5 µg	yellow/gray	✓
93646	Clindamycin	10 µg	dark blue/white	
67237	Erythromycin	15 µg	green	✓
14683	Fleroxacin	5 µg	beige/white	
69531	Gentamycin	30 µg	red/dark green/white	
55941	Imipenem	10 µg	orange/blue	✓
1403	Kanamycin	30 µg	red	✓
19607	Lincomycin	15 µg	blue/dark blue	
5577	Methicillin	10 µg	red/blue/white	
39782	Nalidixic acid	30 µg	orange	✓
56758	Neomycin	30 µg	dark green/white	✓
44543	Norfloxacin	5 µg	pink/gray	
42079	Ofloxacin	5 µg	red/green	✓
54357	Oxacillin	5 µg	tourquize/white	
16244	Penicillin G	6 µg	blue	✓
11435	Pipemidic acid	20 µg	beige/orange	
56383	Piperacillin	100 µg	green/red/blue	✓
56973	Rifampicin	5 µg	red/dark green	✓
75139	Streptomycin	30 µg	pink	
75141	Tetracycline	30 µg	yellow	✓
73477	Trimethoprim	5 µg	gray/white	✓
74794	Trimethoprim + Sulfomethoxazole	1.25 + 23.75 µg	gray/orange/white	
75156	Vancomycin	30 µg	beige/white/dark blue	✓

Table Range of Antimicrobial Susceptibility Disks



Picture
Antimicrobial Susceptibility Disks

Your Sigma-Aldrich Service Partners

Europe

... Austria Wien

Tel. 01-605 81 10
Fax 01-605 81 20
E-mail: sigma@sigma.co.at

... Belgium/Luxembourg Bornem

Free Tel. 0800 147 47
Free Fax 0800 147 45
E-mail: becustsv@europe.sial.com

... Czech Republic Praha

Tel. 246 003 251
Fax 246 003 291
E-mail: CZECustSV@europe.sial.com

... Denmark Brøndby

Tel. 43 56 59 10
Fax 43 56 59 05
E-mail: denorder@europe.sial.com

... Finland Helsinki

Tel. 09-350 9250
Fax 09-350 92555
E-mail: finorder@europe.sial.com

... France St. Quentin Fallavier

Tél. 0800 21 1408 (appel gratuit)
Fax 0800 03 1052 (appel gratuit)
E-mail: fradvsv@europe.sial.com

... Germany Taufkirchen

Free Tel. 0800 51 55 000
Free Fax 0800 64 90 000
E-mail: deorders@europe.sial.com

... Greece Ilioupoli, Athens

Tel. 210-994 8010
Fax 210-994 3831
E-mail: GRCustSV@europe.sial.com

... Hungary Budapest

Tel. (06-1) 269-6474
Fax (06-1) 235-9068
Email: info@sigma.sial.hu

... Ireland Dublin

Tel. 01-404-1900
Fax 01-404-1910
E-mail: EICustsv@europe.sial.com

... Israel Rehovot

Tel. 08-9484-222
Fax 08-9484-200
E-mail: sigisr@sigma.co.il

... Italy Milano

Tel. 02-33417-310
Fax 02-38010-737
E-mail: itorder@europe.sial.com

... Norway Oslo

Tel. 23 17 60 00
Fax 23 17 60 10
E-mail: nororder@sial.com

... Poland Poznań

Tel. 061-829 01 00
Fax 061-829 01 20
E-mail: plcustsv@europe.sial.com

... Portugal Sintra

Tel. 800 20 21 80 (Gratuito)
Fax 800 20 21 78 (Gratuito)
E-mail: poorders@europe.sial.com

... Spain Tres Cantos, Madrid

Tel. 900 10 13 76 (Gratuito)
Fax 900 10 20 28 (Gratuito)
E-mail: esorders@europe.sial.com

... Sweden Stockholm

Tel. 020-35 05 10
Fax 020-35 25 22
E-mail: sweorder@europe.sial.com

... Switzerland Buchs

Free Tel. 0800 80 00 80
Free Fax 0800 80 00 81
E-mail: Fluka@sial.com

... The Netherlands Zwijndrecht

Free Tel. 0800 022 9088
Free Fax 0800 022 9089
Email: nlcustsv@europe.sial.com

... United Kingdom Gillingham

Free Tel. 0800 717181
Free Fax 0800 378785
E-mail: ukcustsv@europe.sial.com

America

... Canada Oakville

Free Tel. 1-800-565-1400
Free Fax 1-800-265-3858
E-mail: canada@sial.com

... USA St. Louis, MO

Toll Free 800 325 3010
Free Fax 800 325 5022
E-mail: custserv@sial.com

Oceania

... Australia Castle Hill

Free Tel. 1800 800 097
Free Fax 1800 800 096
E-mail: ausmail@sial.com

... New Zealand Castle Hill

Free Tel. 0800 936 666
Free Fax 0800 937 777
E-mail: ausmail@sial.com

Africa

... South Africa Johannesburg

Tel. 011-979 1188
Fax 011-979 1119
E-mail: rsa@sial.com

sigma-aldrich.com



The
SIGMA-ALDRICH
Family



SIGMA
Biochemicals and
Reagents for Life
Science Research



ALDRICH
Organics and
Inorganics for
Chemical Synthesis



Fluka
Speciality Chemicals
and Analytical Reagents
for Research



Riedel-de Haën
Laboratory Chemicals
and Reagents for
Research and Analysis



SUPELCO
Chromatography
Products for Analysis
and Purification