

49862 Nitrocefin disks

For the rapid detection of β -lactamase enzymes in isolated colonies of *Neisseria gonorrhoeae*, *Moraxella catarrhalis*, *Staphylococcus spp.*, *Haemophilus influenzae* and anaerobic bacteria.

Composition:

(1 package contains 50 disks in a light resistant plastic vial.)
6mm diameter filter paper disks impregnated with Nitrocefin.

Storage:

Store in the freezer below -10°C in the containers provided. Allow to equilibrate to room temperature before opening then return to freezer storage immediately after use.

Directions:

Place the required number of Nitrocefin disks into a clean empty Petri dish or onto a microscope slide. Disks may be moistened with one drop of deionised water. Do not over-moisten. Using a sterile loop or applicator stick remove several well-isolated and similar colonies and smear onto the surface of a disk. Alternatively: moisten the disk with one drop of deionised water, then holding the disk in forceps, wipe across a colony on an agar plate. Observe the inoculated disk for the development of a red colour.

Interpretation of results:

Positive - Development of a red colour in the area of the disk where the culture was applied. Note the colour change does not normally develop over the whole of the disk.
Negative - No colour change.

A positive result should be interpreted as resistance to penicillin or cephalosporin activity. Susceptibility should be confirmed by standard growth-dependent susceptibility testing methods. Negative results imply but do not guarantee susceptibility.

Quality control:

The list below illustrates a range of performance control strains in routine use:

Test Organisms (ATCC)	Result
<i>Haemophilus influenzae</i> (35036)	positive
<i>Neisseria gonorrhoeae</i> (31426)	positive
<i>Staphylococcus aureus</i> (11632)	positive
<i>Escherichia coli</i> (25922)	negative

User quality control:

Check for signs of deterioration. Quality control must be performed with at least one organism to demonstrate a positive reaction and at least one organism to demonstrate a negative reaction. Do not use the product if the reactions with the control organisms are incorrect.

Limitations:

It is recommended that biochemical and/or serological tests are performed on colonies from pure culture to confirm identification.

For most bacterial strains a positive results will develop within 5 minutes. However, positive reactions for some staphylococci and anaerobic species may take up to 60 minutes to develop.

Detection of staphylococcal β -lactamase is enhanced by testing growth from a medium containing sub-inhibitory concentrations of a β -lactam antibiotic .

