

## 67886 PYrase Strips (Pyrrolidonyl Peptidase Strips)

Diagnostic test for the rapid differentiation of Enterococci from the group D Streptococci and differentiation of *Streptococcus pyogenes* from other haemolytic Streptococci.

### Composition:

(1 package contains 50 test strips and 1.0 ml Reagent)

The kit contains plastic strips with a paper zone saturated with chromogenic substrate for the detection of pyrrolidonyl peptidase and a bottle containing a developing reagent.

### Storage:

Store dry at +2 to +8°C. Expiration can be prolonged when stored at -20°C.

### Directions:

Wipe off several suspect colonies from a Petri dish by the paper zone of diagnostic strip.

Apply approximately 10 µl of developing reagent to the paper zone by plastic loop or micropipette tip after a minute.

The reaction is observed within 1 minute.

### Interpretation of results:


Negative reaction: no colour change develops at the position of wiped colony.

Positive reaction: a red spot develops at the position of wiped colony.

Use Pyrase test due to the recommendation of NRL for Streptococci and Enterococci:

### Result of the test:

Positive (change to red colour)	Negative (no colour change)
Enterococcus	<i>Streptococcus bovis</i>
<i>Streptococcus pyogenes</i>	<i>Streptococcus equinus</i>
<i>Lactococcus lactis</i>	Beta Minute haemolytic Streptococci:
Aerococcus	<i>Streptococcus anginosus</i>
Gamella (most strain)	<i>Streptococcus intermedius</i>
	<i>Streptococcus constellatus</i>
	Other <i>Lactococcus</i> species
	Leuconostoc
	Pediococcus
	Lactobacillus (sporadic +)
	Other species of viridising Streptococci



## References:

1. Godsey et al., Abstr. ASM Annu. Meet., C84, 276 (1981)
2. R.R. Facklam, L.G. Thacker, B. Fox, L. Eriquez, Presumptive identification of streptococci with a new test system, J. Clin. Micro., 15, 987 (1982)
3. B.L. Wasilauskas, K.D. Hampton, Evaluation of the Strep-A-Fluor identification method for strep A streptococci, J. Clin. Micro. 20, 1205 (1984)
4. D.M. Yajko, J. Lawrence, P. Nassos, J. Young, K.W. Hadley, Clinical trial comparing bacitracin with Strp-A-Chek for accuracy and turnaround time in the presumptive identification of *Streptococcus pyogenes*, J. Clin. Micro., 24, 431 (1986)
5. L.P. Gordon, M.A.S. Damm, J.D. Anderson, Rapid presumptive identification of streptococci directly from blood cultures by serologic tests and the L-pyrrolidonyl- $\beta$ -naphthylamide reaction, J. Clin. Micro. 25, 238 (1987)
6. P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, R.H. Tenover, Eds. Manual of clinical microbiology, 7<sup>th</sup> ed., Washington, D.C., ASM Press (1999)
7. J.F. MacFaddin, Biochemical tests for the identification of medical bacteria, 3<sup>rd</sup> ed., Philadelphia, Lippincott Williams & Wilkins (2000)