

1.13301.0001

Microbiology**Bactident® Aminopeptidase****for the detection of L-alanine aminopeptidase in microorganisms****Contents: 50 test strips****Composition**

The reaction zone of a test strip contains:
L-alanine-4-nitroanilide 0,5 µmol; buffering agents

Principle

L-alanine aminopeptidase is an enzyme which is localized in the cell envelope of bacteria and which is found in relevant activities almost exclusively in Gram-negative microorganisms. This enzyme splits off the amino acid L-alanine from various substrates. In the case of these test strips, the substrate L-alanine-4-nitroanilide is split into 4-nitroaniline and the amino acid L-alanine in the presence of alanine aminopeptidase. The presence of the L-alanine aminopeptidase is indicated by the yellow colouration of the 4-nitroaniline. The results of the investigations performed so far indicate that there is a very good correlation between the aminopeptidase reaction and the Gram-behaviour of the microorganisms.

Use

Suspend a thickly grown individual colony (about 2 mm Ø) in 0.2 ml of distilled water to give an opalescent mixture.

Note: Only bacterial colonies without strong intrinsic colourations should be used for the aminopeptidase test.

We recommend that a control test with an aminopeptidase-positive bacterium (e.g. *E.coli*) and an aminopeptidase-negative bacterium (e.g. *Staphylococcus aureus*) should always be carried out at the same time as the main test.

Stability

see expiring date


Only remove the amount of strips needed at the time!
Do not touch the reaction zones of the test strips.
Close receptacle firmly immediately after use. Please store at the specified temperature.

Safe removal


The test strip is to be removed safely after use like bacteria containing material. This may be done by burning, autoclaving or by placing into a 5 to 6 % disinfectant solution – for at least 6 hours.

Procedure


1. Using an inoculating loop, remove an individual, thickly-grown colony from the nutrient medium (fig. 1).




2. In a small test tube, suspend the bacterial mass in 0.2 ml of distilled water (fig. 2).



3. Insert the aminopeptidase test strip into the test tube such that the reaction zone is completely immersed in the bacteria suspension (fig. 3).



4. Incubate the test tube in a water bath (or incubation cupboard) for 10 to a maximum of 30 minutes* at 37 °C.



5. Read off the reaction by comparison with the colour scale (fig. 4).

* A clear yellow colouration of the bacteria suspension can be seen after only 10 minutes in the case of most aminopeptidase-positive microorganisms; if no yellow colouration appears within this time, the incubation should be extended to a maximum of 30 minutes so that the weakly aminopeptidase-positive strains can be recognized or the absence of Gram-negative microorganisms can be confirmed (see table for exceptions).

Aminopeptidase-positive strains*

all Gram-negative microorganisms

The suspension of bacteria turns yellow if L-alanine aminopeptidase-positive organisms are present.

Exceptions: *Bacteroides vulgatus*, *Bacteroides fragilis*, *Campylobacter* species, *Veillonella parvula*

Aminopeptidase-negative strains*

all Gram-positive microorganisms

* acc. to the results of investigations performed so far.

Version 2014-02-05

Merck KGaA, 64271 Darmstadt, Germany
Tel.: +49 (0)61 51 72-24 40

EMD Millipore Corporation, 290 Concord Road
Billerica, MA 01821, USA, Tel. +1-978-715-4321
www.your-contacts.com/biomonitoring

