


# New Products from Aldrich



## Catalysts

Polymer-supported catalysts are easily removed from reaction media.<sup>1</sup>

47,454-1	NAFION® SAC-13 .....	10g; 25g
46,518-6	Palladium, polymer-bound, 10 wt. % Pd .....	1g; 10g

NAFION is a registered trademark of E.I. du Pont de Nemours & Co., Inc.

## Acrylic Monomers

These monomers are used to modify the  $T_g$  and add active or latent functional groups.<sup>2-6</sup>

47,514-9	3-(Trimethoxysilyl)propyl acrylate, 92% .....	5mL; 25mL
47,317-0	[3-(Methacryloylamino)propyl]dimethyl(3-sulfopropyl)ammonium hydroxide salt, 96% .....	50g; 250g
47,096-1	2,2,3,3,3-Pentafluoropropyl acrylate, 98% .....	5mL; 25mL
46,979-3	Poly(propylene glycol) methacrylate, $\bar{M}_n$ ca. 430 .....	100mL; 250mL
46,336-1	2-Butoxyethyl methacrylate, 98% .....	100mL; 250mL
46,106-7	<i>N</i> -(Butoxymethyl)acrylamide .....	100mL

## Amine-Functional Monomers

These materials make thermally stable, high-performance polymers.<sup>7-9</sup>

47,636-6	<i>o</i> -Tolidine sulfone, 99%, mixture of 2,8- and 4,8-dimethyl isomers .....	10g; 50g
47,635-8	4,4'-(1,1'-Biphenyl-4,4'-diyldioxy)dianiline, 98% .....	5g; 50g
47,632-3	4,4'-(1,3-Phenylenedioxy)dianiline, 98% .....	10g; 100g
47,631-5	4,4'-Diaminobenzanilide, 98% .....	10g; 100g
47,630-7	3,4'-Oxydianiline, 97% .....	10g; 50g

## Reactive Polymers

These materials have functionality for crosslinking and chain extension.<sup>7,10,11</sup>

47,457-6	Poly(4-vinylphenol- <i>co</i> -methyl methacrylate) .....	50g; 250g
47,172-0	Oxazoline, polymer-supported, 40 wt. % solution in 1-methoxy-2-propanol/water .....	50mL; 250mL
46,813-4	Poly(styrene- <i>co</i> -allyl alcohol) .....	500g; 1kg
46,833-9	Poly(ethylene- <i>co</i> -1,2-butylene) diol .....	100g; 500g

References: (1) Bucsi, I.; Olah, G.A. *J. Catal.* **1992**, *137*, 12. (2) Moszner, N. et al. *Polym. Bull.* **1996**, *37*, 289. (3) Lee, W.-F.; Tsai, C.-C. *Polymer* **1994**, *35*, 2210. (4) Lacroix Desmazes, P.; Guyott, A. *Polym. Bull. (Berlin)* **1996**, *37*, 183. (5) DeVries, R.A. et al. *Chem. Ind. (Dekker)* **1994**, *53*, 305. (6) Xu, Z. et al. *Zhanjie* **1992**, *13*, 1; *Chem. Abstr.* **1993**, *119*, 204465r. (7) Tamai, S. et al. *Polymer* **1996**, *37*, 3683. (8) Hsiao, S.-H. et al. *J. Polym. Res.* **1995**, *2*, 1. (9) Yang, C.-P. et al. *J. Polym. Res.* **1994**, *1*, 43. (10) Tan, N.C. et al. *Polymer* **1996**, *37*, 3509. (11) Feng, H. et al. *Eur. Polym. J.* **1997**, *33*, 743.



**ALDRICH®** chemists helping chemists in research & industry

P.O. Box 355, Milwaukee, WI 53201 USA Telephone: 414-273-3850 • 800-558-9160 Fax: 414-273-4979 • 800-962-9591 Web Site: [www.sial.com/aldrich](http://www.sial.com/aldrich)

Aldrich is a member of the Sigma-Aldrich family.