


NEW! Aldrich® FT-IR Library on CD-ROM



Quality data
at a low price!

The Aldrich® FT-IR library contains 8,000 condensed and 3,000 vapor phase IR spectra of pure compounds highlighting the most common chemicals found in the Aldrich Handbook. All spectra were acquired, examined, and collected by Sigma-Aldrich, and can be used with Thermo Electron SpecID® (Z40,748-8) spectral search software.

- Nujol and KBR pellets used
- Spectra archived using 2cm^{-1} providing a 4cm^{-1} wave number resolution
- 32-bit ordinate precision
- Frequency range: $400\text{-}4,000\text{cm}^{-1}$

Library includes these chemical categories:

Alkynes	Organometallics
Aromatics	Polymers
Boranes	Steroids
Deuterated compounds	Silanes
Indicators & dyes	And more...
Non-aromatics	

Z54,752-2 \$5,000.00

Aldrich is a registered trademark of Sigma-Aldrich Biotechnology, L.P.
SpecID is a trademark of Thermo Electron Corporation



Aldrich®/SensIR® ATR-IR Library on CD-ROM
The New Look for Fourier Transform Spectra!

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE

SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA
Phone: 800-558-9160 Fax: 800-962-9591 Technical Service: 800-231-8327


SIGMA-ALDRICH