

17776 Chlortetracycline Selective Supplement

An antibiotic supplement recommended for the isolation of *Pencillium viridicatum* and related species isolated from foods.

Composition:
(per vial, sufficient for 1000 ml medium)
Chlortetracycline 50.00 mg

Directions:
Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add to 1000 ml of sterile, molten Pentachloro Rose Bengal Yeast Extract Agar (Fluka No. 17193). Mix well and aseptically pour into sterile petri plates.