

92301 *Listeria mono Selective Supplement I*

A selective supplement recommended by ISO Committee for the isolation of *Listeria* species.

Composition

(per vial, sufficient for 500 ml medium)

Polymyxin B sulphate 38350 IU

Directions

Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water. Mix well and aseptically add it to 460 ml of sterile, molten L. mono Differential Agar Base (Fluka 77408) along with sterile contents of L. mono Enrichment Supplement I (Fluka 03708) and sterile rehydrated contents of L. mono Selective Supplement II (Fluka 91603) or add in 470 ml of sterile, molten L. mono confirmatory Agar Base (Fluka 92302) alongwith sterile content of L. mono Enrichment Supplement II (Fluka 15895) and sterile rehydrated content of L. mono Selective Supplement II (Fluka 91603). Mix well and pour in sterile petri plates.