

Certificate of Analysis

Product-No: 95451

Product: TOC System Suitability Test – Kit

Lot: 1260184 Filling Code 30206921

The Kit contains following products:

86785 System Suitability Solution (1 x 40 ml)
prepared according to the procedure described in USP <643> from lot G1B145 1,4 Benzoquinone

92068 Standard Solution (1 x 40 ml)
prepared according to the procedure described in USP <643> from lot H0B002 Sucrose

89897 Reagent Water (3 x 40 ml)
TOC: < 0.1 mg/l C

Application test according USP <643>: corresponds

Response efficiency: 101.6%

Suitable for all instruments except:

Anatel A-643

Sievers Model 400 (US) or Model 2244 OL (Europe)

The solutions are stabilized and filled thoroughly under laminar flow in pre-cleaned vials.

Storage: The Kit is delivered at 4°C and should be stored in a refrigerator at 4°C until use.

Expiry date: printed on the label

Sigma-Aldrich Production GmbH
Quality Control



Dr. K.D. Schmidt