

Certificate of Analysis

Product-No: 95451
Product: TOC System Suitability Test – Kit

Lot **BCBF5825** **P Code** **11411P01**

The Kit contains following products:

86785 System Suitability Solution (1 x 40 ml)
prepared according to the procedure described in USP <643> from lot H1G332 1,4-Benzoquinone

92068 Standard Solution (1 x 40 ml)
prepared according to the procedure described in USP <643> from lot I0I179 Sucrose

89897 Reagent Water (3 x 40 ml)
TOC: < 0.1 mg/l C

Application test according USP <643>: corresponds

Response efficiency: 96.9% ± 2.2%

Suitable for all instruments except:

Anatel A-643

Sievers Model 400 (US) or Model 2244 OL (Europe)

The solutions are stabilized and filled in pre-cleaned vials.

Storage: The Kit is delivered at 4°C and should be stored in a refrigerator at 4°C until use.

Expiry date: JAN/12

Sigma-Aldrich Production GmbH
Quality Control



Dr. K.D. Schmidt