

アセンティス

エクスプレス

Ascentis[®] Express ガードカラム

2010年4月より発売開始!



- ・カートリッジ充填剤は、2.7 μmのFused-Core™粒子
- ・耐圧9000 psi (62 MPa, 621 bar, 633 kgf/cm²)
- ・手締めで接続可能

Ascentis Express ガード ホルダー

品名	CAT.NO.	価格(¥)
Ascentis Expressガードホルダー(チタニウム製Hybridフェラル付)	53500-U	38,000

Ascentis Express ガード カートリッジ(長さ 5 mm)

固定相	I.D. (mm)	入数	CAT.NO.	価格(¥)
C18	2.1	3	53501-U	50,000
C18	3.0	3	53504-U	50,000
C18	4.6	3	53508-U	50,000
C8	2.1	3	53509-U	50,000
C8	3.0	3	53511-U	50,000
C8	4.6	3	53512-U	50,000
Phenyl-Hexyl	2.1	3	53524-U	50,000
Phenyl-Hexyl	3.0	3	53526-U	50,000
Phenyl-Hexyl	4.6	3	53531-U	50,000
HILIC	2.1	3	53520-U	50,000
HILIC	3.0	3	53521-U	50,000
HILIC	4.6	3	53523-U	50,000

固定相	I.D. (mm)	入数	CAT.NO.	価格(¥)
RP-Amide	2.1	3	53514-U	50,000
RP-Amide	3.0	3	53516-U	50,000
RP-Amide	4.6	3	53519-U	50,000
F5	2.1	3	53594-U	50,000
F5	3.0	3	53597-U	50,000
F5	4.6	3	53599-U	50,000
Peptide ES-C18	2.1	3	53536-U	50,000
Peptide ES-C18	3.0	3	53537-U	50,000
Peptide ES-C18	4.6	3	53542-U	50,000

Ascentis is a registered trademark of Sigma-Aldrich Biotechnology LP. Fused-Core is a trademark of Advanced Materials Technology, Inc. The Optimize Technologies EXP Titanium Hybrid Ferrule and EXP Holders are Patent Pending. The OPTI- prefix is a registered trademark of Optimize Technologies, Inc.

・本カタログに掲載の製品及び情報は2010年4月1日現在の内容であり、掲載の品目、製品情報、価格等は予告なく変更される場合がございます。最新の情報は、弊社Webサイト (sigma-aldrich.com/japan) をご覧ください。
 ・掲載価格は希望納入価格(税別)です。詳細は販売代理店様へご確認ください。
 ・弊社の試薬は試験研究用のみを目的として販売しております。医薬品、家庭用その他試験研究以外の用途には使用できません。

SIGMA-ALDRICH[®]

シグマ アルドリッチ ジャパン株式会社
 アナリティカル事業部

〒140-0002 東京都品川区東品川2-2-24 天王洲セントラルタワー4F

製品に関するお問合せは、弊社テクニカルサポートへ

TEL : 03-5796-7350 FAX : 03-5796-7355

E-mail アドレス : sialjpsp@sial.com

在庫照会/ご注文方法に関するお問合せは、弊社カスタマーサービスへ

TEL : 03-5796-7320 FAX : 03-5796-7325

大阪営業所 : 〒532-0004 大阪市淀川区西宮原2-7-38 新大阪西浦ビル

TEL : 06-6397-5963 FAX : 06-6397-4649

URL : <http://www.sigma-aldrich.com/japan>

お問合せは下記代理店へ

LC 100331 2010.4