

Pure Solv™ Selection Guide

HPLC Grade Solvent	Cat. No.	Cat. No. Pure Solv Micro System		Replacement purification column	Cat. No.
		24/40 Joint	29/32 Joint		
Acetone	270725	Z568031	Z683981	Pure Solv Media (Proprietary)	Z568074
Acetonitrile	34851	Z568031	Z683981	Pure Solv Media (Proprietary)	Z568074
Benzene	270709	Z568023	Z683973	Alumina	Z568066
Chlorobenzene	270644	Z568023	Z683973	Alumina	Z568066
Chloroform	34854	Z568023	Z683973	Alumina	Z568066
Cyclohexane	34855	Z568023	Z683973	Alumina	Z568066
Dichloromethane	34856	Z568023	Z683973	Alumina	Z568066
Diethyl Ether	309966	Z568023	Z683973	Alumina	Z568066
Diisopropylethylamine	496219	Z568031	Z683981	Pure Solv Media (Proprietary)	Z568074
DME	307432	Z568023	Z683973	Alumina	Z568066
DMF	270547	Z568058	Z684007	Molecular Sieve	Z568082
DMSO	34869	Z568023	Z683973	Alumina	Z568066
Ethanol	459828	Z568031	Z683981	Pure Solv Media (Proprietary)	Z568074
Heptane	34873	Z568023	Z683973	Alumina	Z568066
Hexane	293253	Z568023	Z683973	Alumina	Z568066
Methanol	34860	Z568031	Z683981	Pure Solv Media (Proprietary)	Z568074
Pentane	34956	Z568023	Z683973	Alumina	Z568066
Pyridine	270407	Z568031	Z683981	Pure Solv Media (Proprietary)	Z568074
THF	34865	Z568023	Z683973	Alumina	Z568066
Toluene	34866	Z568023	Z683973	Alumina	Z568066
Triethylamine	471283	Z568031	Z683981	Pure Solv Media (Proprietary)	Z568074

SIGMA-ALDRICH
Labware

76385-511653
1051

©2011 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip. Pure Solv is a trademark of Innovative Technology, Inc.

SIGMA
Life Science

ALDRICH
Chemistry

Fluka
Analytical

SUPELCO
Analytical

SIGMA-ALDRICH®