



SIGMA-ALDRICH

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

N-ACETYL-L-CYSTEINE

Sigma Prod. No. **A7250**,
Fluka 01039, Aldrich 138061

CAS NUMBER: 616-91-1

SYNONYMS: L-alpha-acetamido-beta-mercaptopropionic acid, acetin, acetylcysteine, N-acetylcysteine, N-acetyl-N-cysteine, N-acetyl-3-mercaptoalanine, airbron, broncholsin, flumucetin, flumucil, flumicil, inspir, mercapturic acid, mucolyticum, mucolyticum lappe, mucolytikum lappe, mucomyst, mucosolvin, NAC, NAC-TB, NSC 111180, parvolex, respaire

PHYSICAL PROPERTIES:

Appearance: White to white with light yellow cast powder

Melting Point: 109-110°C¹

Molecular formula: C₅H₉NO₃S

Formula weight: 163.2 (anhydrous)

pKa: 9.5 at 30°C²

Optical rotation: +5°(c = 3% in water)³

Purity: Not less than 99% (TLC)

METHOD OF PREPARATION:


Sigma's product is prepared from the acetylation of cysteine from human hair.

STABILITY / STORAGE AS SUPPLIED:

This product if stored at 2-8°C has a shelf-life of 3 years.

SOLUBILITY / SOLUTION STABILITY:

N-acetyl-L-cysteine is soluble 1 in 8 of water and 1 in 2 of ethanol. It is practically insoluble in chloroform and ether.² Sigma has dissolved this product in water at 100 mg/mL with heating and obtained a clear, colorless solution.



N-ACETYL-L-CYSTEINE
Sigma Prod. No. **A7250**,
Fluka 01039, Aldrich 138061

SOLUBILITY / SOLUTION STABILITY: (continued)

Aqueous solutions of cysteine oxidize to cystine on contact with air at neutral or alkaline pH. Oxidation is accelerated by traces of heavy metals, especially copper and iron; relatively stable in acid.⁴ This information is for Cysteine but will likely be applicable to N-acetyl-L-cysteine as well.

APPLICATIONS:

Some applications of N-acetyl-L-cysteine include its use as an anti-inflammatory drug to prevent liver damage by tylenol overdoses,⁵ as a sputum liquefying agent,^{6,7} and as an inhibitor of HIV replication.⁸ A review of animal pharmacology and biochemistry is presented in *Pharmacological and Biochemical Properties of Drug Substances*, vol. 2, 479-488 (M.E. Goldberg, Ed., American Pharm. Association, Washington, D.C., 1979).

REFERENCES:

1. *The Merck Index*, 11th ed., p. 14, #82 (1989).
2. *Clarke's Isolation and Identification of Drugs*, 2nd ed., p. 318 (1986).
3. *Dictionary of Organic Compounds*, 5th ed., vol. 1, p. 35 (1982).
4. *Data for Biochemical Research*, 3rd ed., p. 12-13 (1986).
5. Lauterburg, B.H., et al., *J. Clin. Invest.*, **71**, 980-991 (1993).
6. Kubica, G.P., et al., *Amer. Respir. Dis.*, **87**, 775-779 (1963).
7. Hirsch, S.R., et al., *J. Lab. Clin. Med.*, **74**, 346-353 (1969).
8. Kalebic, T., et al., *Proc. Nat. Acad. Sci. USA*, **88**, 986-990 (1991).

Sigma warrants that its products conform to the information contained in this and other Sigma!Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice of packing slip.