


### p3XFLAG-CMV-9 (6.4 kb)



### p3XFLAG-CMV-10 (6.4 kb)



## Multiple Cloning Site

(p3XFLAG-CMV-9\* and p3XFLAG-CMV-10)

### 3XFLAG Peptide Sequence

Met\* Asp Tyr Lys Asp His Asp Gly Asp Tyr Lys Asp His Asp Ile

ATG GAC TAC AAA GAC CAT GAC GGT GAT TAT AAA GAT CAT GAC ATC

TAC CTG ATG TTT CTG GTA CTG CCA CTA ATA TTT CTA GTA CTG TAG

### 3XFLAG Peptide Sequence

Asp Tyr Lys Asp Asp Asp Asp Lys

GAT TAC AAG GAT GAC GAT GAC AAG CTT GCG GCC GCG AAT TCA TCG ATA

CTA ATG TTC CTA CTG CTA CTG TTC GAA CGC CGG CGC TTA AGT AGC TAT

Hind III

Bgl II EcoR V Kpn I Xba I Bam HI

GAT CTG ATA TCG GTA CCA GTC GAC TCT AGA GGA TCC CGG GTG

CTA GAC TAT AGC CAT GGT CAG CTG AGA TCT CCT AGG CCC CAC

\*For pFLAG-CMV-9, the Met-preprotrypsin leader sequence (PPT LS) precedes the FLAG coding sequence.