

Product Information

ANTI-CANNABINOID RECEPTOR 2, HUMAN

Developed in Rabbit, Affinity Isolated Antibody

Product Number **C 1358**

Product Description

Anti-Cannabinoid Receptor 2, Human (CB2) is developed in rabbit using a highly purified fusion protein containing the first 33 amino acid residues of the human CB2 receptor as the immunogen. The antibody is affinity isolated on immobilized immunogen.

Anti-Cannabinoid Receptor 2, Human specifically recognizes the CB2 receptor from human and rat by immunoblotting and immunocytochemical staining of living or fixed cultured hippocampal neurons.

Cannabinoids exert their well known physiological effects through two G-protein coupled receptors, CB1 and CB2.¹ Both receptors have been shown to inhibit adenylyl cyclase as well as stimulate the mitogen-activated protein kinase, MAPK. CB1 receptors also modulate ion channels through direct G-protein interactions. Δ^9 -Tetrahydrocannabinol and related ligands exert their psychoactive effects by inhibiting presynaptic N- and Q-type calcium channels.²

The CB1 receptor is primarily expressed in brain (particularly cortex, striatum, hippocampus and cerebellum), spinal cord, eye and testis.³ CB1 receptors are involved in the modulation of nociception (pain), movement, learning and memory, emotion and many other neuronal processes.⁴

The CB2 receptor is 44 % identical to CB1 at the amino acid level. It is thought to function primarily in the peripheral immune system, having been found in the spleen. It may be present in the central nervous system, including the retina.⁵

Reagent

Anti-Cannabinoid Receptor 2, Human is supplied as 100 μ l of affinity isolated rabbit antibody in phosphate buffered saline containing 1 mg/ml bovine serum albumin, 50 % glycerol and 0.05 % sodium azide.

Precautions and Disclaimer

Due to the sodium azide content, a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling practices.

Storage/Stability

Store the antibody at -20°C . Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

The recommended working dilution is 1:250 for immunoblotting. This antibody recognizes a single band at approximately 60 kDa.

The recommended working dilution is 1:500 for immunocytochemistry using standard techniques.

Note: In order to obtain best results and assay sensitivities of different techniques and preparations, we recommend determining optimal working dilutions by titration test.

References

1. Felder, C. C. and Glass, M., *Ann. Rev. Pharmacol. Toxicol.*, **38**, 179-200 (1998).
2. Twitchell, W., et al., *J. Neurophysiol.*, **78**, 43-50 (1997).
3. Tsou, K., et al., *Neuroscience*, **83**, 393-411 (1998).
4. Pertwee, R.G., *Curr. Med. Chem.*, **6**, 635-664 (1999).
5. Lu, Q., et al., *Vis. Neurosci.*, **17**, 91-95 (2000).

mje 10/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.