

## Binding of Protein L to Various Immunoglobulin Light Chains

Species	Protein L Binding
Human kappa I	++++
Human kappa II	-
Human kappa III	++++
Human kappa IV	++++
Human kappa I-IV	-
Human kappa V	+
Human kappa I	++++
Human kappa II	-
Human kappa V	+

This table is not meant to be complete, but only covers the known kappa subgroups that protein L-binding has been tested on. It should be noted that protein L is restricted to specific subclasses of the VL domain. Thus, the affinity indicated in the table is not general for the IgG subclass, but accounts only for those antibodies carrying the right kappa light chain.