

Product Information

HEMATOXYLIN & EOSIN

(Progressive Method-Gill's Hematoxylin Template)

LEICA ST4040

Standard Version

Procedure GHS

Begin	Xylene	Xylene	Xylene	Xylene	100% Alcohol	100% Alcohol	100% Alcohol	95% Alcohol	Wash	Wash	Gill's Hem 2	Gill's Hem 2	Gill's Hem 2	Wash
Wash	Wash	Scott's Tap	Wash	Wash	Eosin Y	95% Alcohol	100% Alcohol	100% Alcohol	100% Alcohol	Xylene	Xylene	Xylene	End	