

## Product Information

### PAPANICOLAOU (PAP) STAIN FOR CYTOLOGIC SPECIMEN

(Progressive Method-Gill's Papanicolaou Template)

### LEICA ST4040

Procedure PAP-GHS

Begin	95% Alcohol	80% Alcohol	Wash	Wash	Gill's Hem 2	Gill's Hem 2	Gill's Hem 2	Wash	Wash	Scott's Tap	Wash	Wash	OG-6	OG-6
OG-6	OG-6	100% Alcohol	EA-50	EA-50	EA-50	EA-50	EA-50	100% Alcohol	100% Alcohol	100% Alcohol	Xylene	Xylene	End	