

## Product Information

### HEMATOXYLIN & EOSIN

(Progressive Method-Gill's Hematoxylin Template)

### LEICA AutoStainer XL Stainer

Procedure GHS

Station 1	Station 2	Station 3	Station 4	Station 5	Station 6	Station 7	W1	W2	W3	W4	W5	
<b>Xylene</b> 1min, 30 sec	<b>Xylene</b> 1 min, 30 sec	<b>Xylene</b> 1 min, 30 sec	<b>100% Alc</b> 1 min, 30 sec	<b>100% Alc</b> 1 min, 30 sec	<b>95% Alc</b> 30 sec	<b>80% Alc</b> 30 sec	<b>Wash</b> 1 min	<b>Wash</b> 2 min	<b>Wash</b> 2 min, 30 sec	<b>Wash</b> 1min, 30 sec	<b>Wash</b> 1min	<b>Oven</b>
<b>End</b>	Station 18 <b>Xylene</b> 2 min	Station 17 <b>Xylene</b> 2 min	Station 16 <b>100% Alc</b> 1 min, 30 sec	Station 15 <b>100% Alc</b> 1 min	Station 14 <b>100% Alc</b> 30 sec	Station 13 <b>Not Used</b>	Station 12 <b>95% Alc</b> 20 sec	Station 11 <b>Eosin Y</b> 30 sec	Station 10 <b>Not Used</b>	Station 9 <b>Scott's Tap</b> 1 min	Station 8 <b>Gill's 2 Hem</b> 3 min, 30 sec	<b>Hold</b>