


#501413 Ma Quick Connect
w/Vinyl Cap

#500123 Tubing Plug

Double Poly Bagged
w/501754 JRH Label &
#501403 Irradiation Dot

2.5' #501947 1/4" x 3/8" Clear
C-Flex Tubing w/Pinch Clamp

3' #501610 1/8" x 1/4" Clear
C-Flex Tubing w/Pinch Clamp



Each Item Will Be Individually
Double Poly Bagged & All Tubing Will
Be Coiled

#502103 100L JRH Bulk Bag

P1

P2



CONFIDENTIAL; THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH
ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

BIOEAZE™ 100L PE ST BAG TER

Item #1947B