



Microbiology Media & Components

Save 33% on these Microbiology Media & Components:

Cat. No.	Product Description	Price (€)
S-Gal Products		
S-Gal		
S9811-100MG	Powder, replacement for X-Gal; for black color development to occur, ferric ion must be present	18.15
S-Gal sodium salt		
S7313-100MG	Powder, more water soluble version for S-Gal, soluble >50 mg/mL	17.16
S-Gal/LB Agar Blend, without IPTG		
S9938-500ML		34.19
S-Gal/LB Agar Blend		
C4478-500ML	Autoclavable, microwavable growth medium complete with S-Gal and IPTG to make 500 mL	28.25

Cat. No.	Product Description	Price (€)
LB Agar & Broth		
LB Agar		
L2897-1KG	Powder, Tryptone (pancreatic digest of casein)	62.50
L2897-250G	10 g/L, Yeast extract 5 g/L, NaCl 5 g/L, Agar 15 g/L	20.66
LB Agar		
L7025-100TAB	Tablet, 1.68 g per tablet, dissolve the tablet in 50 mL of deionized water.	35.71
L7025-500TAB		134.97
LB Broth		
L3022-1KG	Powder, Tryptone (pancreatic digest of casein)	62.50
L3022-250G	10 g/L, Yeast extract 5 g/L, NaCl 5 g/L	21.32
LB Broth		
L7275-100TAB	Tablet, 1.1 g per tablet. Dissolve the tablet in 50 mL of water.	36.30
L7275-500TAB		124.74
LB Agar EZMix Powder		
L7533-1KG	Granulated, dust free powder, Tryptone (pancreatic digest of casein) 10 g/L,	183.15
L7533-6X500ML	Yeast extract 5 g/L, NaCl 5 g/L, Agar 15 g/L	30.43
LB Broth EZMix Powder		
L7658-1KG	Granulated, dust free powder, Tryptone (pancreatic digest of casein) 10 g/L,	74.25
L7658-6X500ML	Yeast extract 5 g/L, NaCl 5 g/L	17.29

To place an order, please contact your local sales office and quote **promo code X65**.

Offer expires 31 December 2009.