

Quality Starts with Corning® Stripettes®

CORNING

Buy 10 cases of Corning Stripette Serological Pipets and get a free Stripettor™ Pipetting Aid.

(details on back).

Are you spending a lot of time and using a lot of materials re-doing procedures with pipets that are leaking at the seam, dripping from the tip, contaminating your workspace or transferring inaccurate volumes? Then switch to Corning Stripettes and spend less time, waste fewer materials, and have more confidence in your results.

Save time and money by using Corning's high quality Stripettes:

- ▶ RNase-/DNase-free to prevent contamination
- ▶ Nonpyrogenic to promote effective cell growth
- ▶ Color-coded magnifier stripe to easily see meniscus volume
- ▶ Anti-drip tip and strong seal to prevent leaks and assure accurate transfer volumes and uncontaminated workspaces




- RNase-/DNase-free
- Nonpyrogenic



Anti-drip tip




Strong seal to prevent leaks at the seam



Color-coded stripe for magnified volume and meniscus reading




Also available in triple-bagged packaging



SIGMA-ALDRICH®

Buy 10 cases of Corning® Stripette® Serological Pipets and get a free Stripettor™ Pipetting Aid.

Corning Stripette Serological Pipets Ordering Information



Stripette Pipets, Individually Wrapped, Clear Plastic Wrap

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color Coded Stripe	Qty/ Pk	Qty/ Cs
CLS4011	1	1/100	0.2	Yellow	100/bag	1000
CLS4012	1	1/100	0.2	Yellow	100/bag	200
CLS4021	2	1/100	0.2	Green	100/bag	1000
CLS4051	5	1/10	2.5	Blue	50/bag	200
CLS4101	10	1/10	3.0	Orange	50/bag	200
CLS4492*	10	1/10	3.0	Orange	50/bag	200
CLS4251	25	2/10	10.0	Red	50/bag	200
CLS4501	50	1/2	10.0	Purple	25/bag	100
CLS4484	100	1	N/A	Aqua	10/bag	100

*Cat. No. 4492 features a wider tip opening for handling viscous fluids.

Stripette Pipets, Individually Wrapped, Paper/Plastic Wrap

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color Coded Stripe	Qty/ Pk	Qty/ Cs
CLS4485	1	1/100	0.2	Yellow	50/bag	1000
CLS4486	2	1/100	0.2	Green	50/bag	1000
CLS4487	5	1/10	2.5	Blue	50/bag	200
CLS4488	10	1/10	3.0	Orange	50/bag	200
CLS4489	25	2/10	10.0	Red	25/bag	200
CLS4490	50	1/2	10.0	Purple	25/bag	100
CLS4491	100	1	N/A	Aqua	10/bag	100

Stripette Pipets, Bulk Packed in Bags

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color Coded Stripe	Qty/ Pk	Qty/ Cs
CLS4010	1	1/100	0.2	Yellow	50/bag	1000
CLS4020	2	1/100	0.2	Green	50/bag	1000
CLS4050	5	1/10	2.5	Blue	50/bag	500
CLS4100	10	1/10	3.0	Orange	50/bag	500
CLS4250	25	2/10	10.0	Red	25/bag	200
CLS4500	50	1/2	10.0	Purple	25/bag	100

Stripette Pipets, Clean Room Pack, Individually Wrapped, Paper/Plastic, Triple Bagged

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color Coded Stripe	Qty/ Pk	Qty/ Cs
CLS7041	1	1/100	0.2	Yellow	50/bag	1000
CLS7042	2	1/100	0.2	Green	50/bag	1000
CLS7045	5	1/10	2.5	Blue	50/bag	200
CLS7015	10	1/10	3.0	Orange	50/bag	200
CLS7016	25	2/10	10.0	Red	25/bag	200
CLS7017	50	1/2	10.0	Purple	25/bag	100
CLS7000	100	1	N/A	Aqua	10/bag	100

Stripettor™ Pipetting Aid

Cat. No.	Product Description	Qty/Cs
CLS4910	Stripettor with sterile filter, rechargeable battery and recharger/adaptor	1

Redeem the Stripettor Pipetting Aid:

Please quote promotion code W51 when placing your order.

Please do not enter the free stripettor in your order, it will be shipped separately within 1 month.

Contact details: Visit www.sigmaaldrich.com for your local office phone numbers and email addresses.

SIGMA-ALDRICH®

Offer valid until December 31, 2009. Offer void where prohibited by law or company policy. Offer cannot be combined with other Corning Stripette offers.

Corning and Stripette are registered trademarks of Corning Incorporated, Corning, NY. Stripettor is a trademark of Corning Incorporated, Corning, NY. Corning Incorporated, One Riverfront Plaza, Corning, NY 14831-0001

CORNING

Corning Incorporated
Life Sciences

For product information:
t +31 20 659 6051
f +31 20 659 7673

www.corning.com/lifesciences