

## 04327 Atto 612Q

### Application

Atto 612Q is a new quencher for fluorescence characterized by high molecular absorption (114.000). Atto 612Q is characterized by a high photostability and thermostability. It is moderately hydrophilic, and soluble in polar solvents as DMF or DMSO.

### Product Description

$\epsilon$	114.000
$\eta_{fl}$	45% (in water)
Abs. max.:	615nm(H <sub>2</sub> O)
MW:	791
Quantity	1 mg

Absorbance (arbitrary units)

