

# SPE Promotion



Solid phase extraction is a form of digital (step-wise) chromatography designed to extract, partition, and/or adsorb one or more components from a liquid phase (sample) onto stationary phase (sorbent or resin).

Over the last twenty years, SPE has become the most powerful technique available for rapid and selective sample preparation (prep) prior to analytical chromatography. SPE extends a chromatographic system's lifetime, improves qualitative and quantitative analysis, and by changing an analyte of interest's original matrix environment to a simpler matrix more suitable for subsequent analysis, the demand placed on an analytical instrument is considerably lessened.

**For a special time period, we are offering following products on a very special price.**

**If you buy the Manifold and the Vacuum pump, You will get a pack of the SPE tubes of your choice **FREE****

Cat. No.	Product Name	List Price (€)	Special Price (€)
57044	Visiprep™ SPE Vacuum Manifold, DL (Disposable Liner), 12-port model	<del>1015</del>	711
57160-U	Preppy™ 12-Port Vacuum Manifold	<del>535</del>	375
Z288268EU-1EA	KNF Laboport® mini-pump Pump with vacuum gauge and regulator, AC230 V	<del>392</del>	274

## SPE Tubes

52603-U	Discovery® DSC-18 SPE Tube, bed wt. 500 mg, volume 3 mL, pk of 54	<del>134</del>	94
52606-U	Discovery® DSC-18 SPE Tube, bed wt. 1 g, volume 6 mL, pk of 30	<del>100</del>	70
57063	Supelclean™ ENVI-18 SPE Tube, bed wt. 500 mg, volume 3 mL, pkg of 54 ea	<del>133</del>	93
57064	Supelclean™ ENVI-18 SPE Tube, bed wt. 500 mg, volume 6 mL, pkg of 30 ea	<del>70</del>	49
57137	Supelclean™ ENVI-18 SPE Tube, bed wt. 5 g, volume 20 mL, pk of 20	<del>200</del>	196
505706	Supelclean™ ENVI-18 SPE Tube, bed wt. 1 g, volume 6 mL, pk of 30	<del>90</del>	69
57088	Supelclean™ ENVI-Carb™ SPE Tube, bed wt. 250 mg, volume 3 mL, pk of 54	<del>148</del>	104
57094	Supelclean™ ENVI-Carb™ SPE Tube, bed wt. 500 mg, volume 6 mL, pk of 30	<del>136</del>	95
57012	Supelclean™ LC-18 SPE Tube, bed wt. 500 mg, volume 3 mL, pkg of 54 ea	<del>133</del>	93
505471	Supelclean™ LC-18 SPE Tube, bed wt. 1 g, volume 6 mL, pk of 30	<del>90</del>	69
57017	Supelclean™ LC-SAX SPE Tube, bed wt. 500 mg, volume 3 mL, pkg of 54 ea	<del>177</del>	124
57057	Supelclean™ LC-Florisil® SPE Tube, bed wt. 1 g, volume 6 mL, pkg of 30 ea	<del>121</del>	85
57116	Supelclean™ LC-Si SPE Tube, bed wt. 2 g, volume 12 mL, pk of 20	<del>99</del>	69
505048	Supelclean™ LC-Si SPE Tube, bed wt. 500 mg, volume 3 mL, pk of 54	<del>107</del>	75
505374	Supelclean™ LC-Si SPE Tube, bed wt. 500 mg, volume 6 mL, pk of 30	<del>109</del>	76

**Promo code: U77**

Offer valid until: 31.03.2008

**Sigma-Aldrich (O.M.) Ltd.**

Αργοναυτών 72  
163 46 Ηλιούπολη

Τηλ.: 210 99 48 010  
Fax: 210 99 43 831

E-mail: [GRCustSV@europe.sial.com](mailto:GRCustSV@europe.sial.com)  
[www.sigma-aldrich.com](http://www.sigma-aldrich.com)