

Bulb-to-Bulb Distillation in a Flash!

Quickly Distill the Most
Difficult Materials

Adaptable, Safe...Precise



Aldrich® Kugelrohr Distillation Apparatus

Set includes: air-bath oven with digital temperature controller, rotary drive, and glassware set containing a straight tube, 25- and 100-mL oven flasks, 25- and 100-mL single bulb tubes all with 14/20 joints. CE compliant.



Special Offer!

Purchase a Kugelrohr Distillation Apparatus before June 30, 2006 and receive a **FREE** 1-Liter Kugelrohr glassware set with greaseless, System 45® threaded connections. Specify Catalog No. **Z510181** and the Promo Code when placing your order for Kugelrohr.

| Voltage | Cat. No. | Promo Code | Price |
|---------|---------------------|------------|------------|
| 115 | Z401137-1SET | PRO100230 | \$3,154.60 |
| 230 | Z401145-1SET | PRO100229 | 3,653.00 |

New
Glass Center



sigma-aldrich.com/glass

See it...Select it...Have it
tomorrow! It's that simple.

For more information, please visit sigma-aldrich.com/labequipment.

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA

Aldrich and Precision-Seal are registered trademarks of Sigma-Aldrich Biotechnology, L.P. System 45 is a registered trademark of NDS Technologies, Inc.



SIGMA-ALDRICH