



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

**ANTI-HUMAN IgM (m-CHAIN SPECIFIC)
PEROXIDASE CONJUGATE**
Antibody developed in Goat
Affinity Isolated Antigen Specific Antibody

Product Number **A 0420**

Product Description

Anti-human IgM (μ -chain specific) is developed in goat using the human IgM as the immunogen. The antibody is isolated from goat anti-human IgM antiserum by immunospecific purification to remove essentially all goat serum proteins, including immunoglobulins, which do not specifically bind to the μ -chain of human IgM. Goat anti-human IgM is conjugated to peroxidase by means of a two-step glutaraldehyde method.

Specificity of the Peroxidase Conjugated Anti-Human IgM is determined by Enzyme Linked Immunosorbent Assay (ELISA). The conjugate is specific for human IgM when tested against human IgA, IgG, IgM, Bence Jones kappa and lambda myeloma proteins.

Cross-reactivity of the antibody-conjugate is determined by ELISA. The conjugate shows no reactivity with mouse or rat IgG.

Identity and purity of the antibody is established by immunoelectrophoresis (IEP), prior to conjugation. Electrophoresis of the antibody preparation followed by diffusion versus anti-goat IgG and anti-goat whole serum results in single arcs of precipitation. The product is purified to remove unconjugated material.

Reagents

The product supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 0.01% thimerosal as a preservative.

Precautions and Disclaimer

Consult the MSDS for information regarding hazards and safe handling practices.

Storage/Stability

Store at 2-8 °C for up to one month.

For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Product Profile

Molar Ratio (IgG:Peroxidase) = 0.8 to 1.5

Titers

1. Direct ELISA: Minimum 1:50,000
We are now reporting lot specific information as a titer by direct ELISA rather than as a working dilution (see below). Titer is defined as the dilution of conjugate sufficient to give a change in absorbance of 1.0 at 450 nm after 30 minutes of substrate conversion at 25 °C (Voller, et al.)¹. Microtiter plates are coated with purified human IgM at a concentration of 5 μ g/ml in 0.05 M carbonate-bicarbonate buffer, pH 9.6 (Carbonate-Bicarbonate Bicarbonate Buffer Capsules are available as Sigma Product No. C 3041).

Substrate: o-Phenylenediamine Dihydrochloride (OPD, Sigma Product No. P 8287), 0.4 mg/ml in 0.05 M phosphate-citrate buffer, pH 5.0 containing 0.03% sodium perborate (Phosphate-Citrate Buffer Capsules with Sodium Perborate are available as Sigma Product No. P 4922).

2. Dot Blot
 - a. A minimum dilution of 1:6,000 was determined by direct dot blot assay using 40 ng human IgM/dot.
 - b. A minimum dilution of 1:100,000 was determined in a direct chemiluminescence assay using 20 ng human IgM/dot. Luminol plus enhancer was used as substrate.
3. Immunohistology
A minimum dilution of 1:200 was determined in a direct assay using formalin-fixed, paraffin-embedded human tonsil sections.

Working dilutions should be determined by titration assay. Due to differences in assay systems, these titers may not reflect the user's actual working dilution.

This goat antisera was maintained at pH 5.0 for 40 minutes to meet USDA requirements.

References

1. Voller, A., et al., Bull. World Health Organ., **53**, 55 (1976).

JWM/DAA 12/2002

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.