



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## Product Information

### Monoclonal Anti-CREB

#### Clone 48H2

produced in rabbit, purified immunoglobulin

Catalog Number **C8977**

### Product Description

Monoclonal Anti-CREB (rabbit IgG isotype) is derived from the 48H2 hybridoma produced by immunizing rabbits with GST-CREB full length fusion protein. The antibody is purified using Protein G chromatography.

Monoclonal Anti-CREB detects endogenous levels of total CREB protein. The antibody reacts with human, mouse, rat and monkey CREB. It may be used for immunoblotting, immunoprecipitation, immunocytochemistry, immunofluorescence, and immunohistochemistry.

CREB is a bZIP transcription factor that binds the cyclic AMP response element (CRE) and activates transcription in response to a variety of extracellular signals including cAMP, membrane depolarization, increased intracellular  $\text{Ca}^{2+}$ , and growth and neurotrophic factors.<sup>1-3</sup> Phosphorylation of CREB at Ser<sup>133</sup> regulates the ability of CREB to activate transcription when bound at a CRE. Mutation of Ser<sup>133</sup> renders CREB non-responsive to multiple signaling pathways. A variety of protein kinases have been shown to phosphorylate CREB at Ser<sup>133</sup> *in vitro*, including PKA, PKC, CAM kinase II and IV, and p90<sup>sk</sup>. Phosphorylation at Ser<sup>133</sup> has also been shown to promote interaction with a CREB binding protein, CBP, required for transcriptional activation by CREB, AP-1, and SRF dependent promoters.<sup>3,4</sup> CREB appears to play an important role in learning and memory in both flies and mice. Mice lacking CREB exhibit deficiencies in spatial learning tasks, while flies overexpressing or lacking CREB show enhanced or diminished learning.<sup>5,6</sup>

### Reagent

Supplied as a solution in 10 mM HEPES sodium, pH 7.5, containing 150 mM sodium chloride, 100 µg/mL bovine serum albumin, 50% glycerol and <0.02% sodium azide.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

### Storage/Stability

Store at  $-20^{\circ}\text{C}$ . Do not aliquot the antibody. Working dilution samples should be discarded if not used within 12 hours.

### Product Profile

Immunoblotting: recommended working dilution is 1:1,000 using an extract from SK-N-MC, NIH-3T3, and COS cells.

Immunoprecipitation: recommended working dilution is 1:250.

Immunocytochemistry (indirect immunofluorescence): recommended working dilution is 1:1,600.

Immunohistochemistry: recommended working dilution is 1:400 using paraffin-embedded sections or frozen sections.

Flow cytometry: recommended working dilution is 1:400 using SK-N-MC cells.

**Note:** In order to obtain the best results in various techniques and preparations, we recommend determining optimal working dilution by titration.

### References

1. Gonzalez, G.A., and Montminy, M.R., *Cell*, **59**, 675-680 (1989).
2. Sheng, M., et al., *Science*, **252**, 1427-1430 (1994).
3. Kwok, R.P.S., et al., *Nature*, **370**, 223-226 (1994).

4. Arias, J., et al., *Nature*, **370**, 226-229 (1994).
5. Frank, D.A., and Greenberg, M.E., *Cell*, **79**, 5-8 (1994).
6. Yin, J.C.P., et al., *Cell*, **81**, 107-115 (1995).

Rabbit monoclonals produced using Epitomics® Technology,  
U.S. Patent No. 5,675,063.  
Epitomics is a registered trademark of Epitomics, Inc.

KAA,PHC 02/07-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.