



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone (800) 325-5832 (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## Product Information

### ANTI-EPHB1 (ELK RECEPTOR)-SAM

Developed in Sheep, IgG Fraction of Antiserum

Product Number **E 6392**

#### Product Description

Anti-EphB1 (Elk Receptor)-SAM is developed in sheep using a highly purified fusion protein of GST and EphB1 (Elk Receptor)-SAM corresponding to amino acid residues 883-984 of human EphB1 as the immunogen.

Anti-EphB1 (Elk Receptor)-SAM recognizes the sterile alpha motif region of human EphB1 protein by immunoblotting and immunoprecipitation.

EphB family proteins bind ephrins of the B class. EphB1 is expressed predominately in developing neural structures in embryos, and in vascular epithelium of kidney, and other tissues. Upon binding to alternatively oligomerized ephrin B1, EphB1 signals regulation of cell attachment and cell-cell assembly. Members of this protein family are implicated in neuronal and vascular cell targeting.

#### Reagents

Anti-EphB1 (Elk Receptor)-SAM is supplied as 250 µg of purified IgG at a concentration of 0.2 µM in phosphate buffered saline with 0.08% sodium azide.

#### Precautions and Disclaimer

Due to the sodium azide content, a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling practices.

#### Storage/Stability

Antibodies should be stored at -20°C. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

#### Product Profile

The recommended working dilution is 1-10 µg/ml for immunoblotting using peroxidase conjugated anti-sheep IgG and detection by chemiluminescence. The recommended working dilution is 10-20 µg/ml for immunoprecipitation.

Note: In order to obtain best results and assay sensitivities of different techniques and preparations, we recommend determining optimal working dilutions by titration test.

#### References

1. Gao, P.P. et al., J. Neurosci. Res., **60(4)**, 427-436 (2000).
2. Wahl, S. et al., J. Cell Biol., **149(2)**, 263-270 (2000).

mje 8/00

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.