

Product Information

MONOCLONAL ANTI- α -TUBULIN CONJUGATE CLONE DM1A FITC CONJUGATE Purified Mouse Immunoglobulin

Product No. **F 2168**

Product Description

Monoclonal Anti- α -Tubulin (mouse IgG1 isotype) is derived from the DM1A hybridoma produced by the fusion of mouse myeloma cells and splenocytes from immunized BALB/c mice. Purified chick brain microtubules were used as immunogen.¹ The isotype is determined using Sigma ImmunoType™ Kit (Sigma ISO-1) and by a double diffusion immunoassay using Mouse Monoclonal Antibody Isotyping Reagents (Sigma ISO-2). The product is Protein A-purified Monoclonal Anti- α -Tubulin antibody conjugated to fluorescein isothiocyanate, isomer I. It is purified by gel filtration and contains no detectable free FITC.

FITC Conjugated Monoclonal Anti- α -Tubulin specifically recognizes an epitope in the carboxy-terminal part of α -tubulin.² It localizes α -tubulin in human, monkey, bovine, chicken, goat, murine, rat, gerbil, hamster, rat kangaroo, amphibia, sea urchin, trypanosome, yeast, fungi and tobacco.

Monoclonal Anti- α -Tubulin may be used for the detection and localization of α -tubulin using various immunocytochemical and immunohistochemical assays. It may also be used in dual staining assays.

α/β -Tubulin is the major building block of microtubules. These intracellular, hollow, cylindrical, filamentous structures are present in virtually all eukaryotic cells. Self-assembly of α/β -tubulin leads to polar microtubular structures built from linearly arranged strings of alternatively α - and β -tubulin pairs pointing in the same direction. Microtubules function as structural and mobile elements in mitosis, intracellular transport, ciliary flagellar motility and generation and maintenance of cell shape. α/β -Tubulin and γ -tubulin are members of the tubulin superfamily of proteins. α/β -Tubulin is a heterodimer which consists of one α -tubulin chain and one β -tubulin chain; each subunit has a molecular weight of 55 kD and they share considerable homology.^{3,4,5} Both subunits exist in multiple isotypes in many organisms and tissues.⁶

The C-terminal portions of the two subunits have an unusually high proportion of acidic residues. Tubulin can be postranslationally modified by phosphorylation (β -tubulin) or by acetylation, deacetylation and glutamylation (α -tubulin). Tubulin is a highly conserved protein although exceptional cases have been reported. It specifically interacts with ions, nucleotides, microtubule-associated proteins (MAPs) and drugs. The large N-terminal domain of the α and β -tubulin subunits bind GTP. GTP binding to α -tubulin is non-exchangeable while that to the β -subunit is exchangeable. This GTP is hydrolysed upon microtubule assembly. The smaller C-terminal domain of the two subunits bind MAPs. Both domains are involved in microtubular self-assembly. Antibodies to tubulin are specific and useful tools in studying the intracellular distribution of tubulin and the static and dynamic aspects of the cytoskeleton.

Reagents

The product is provided as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 1% BSA and 15mM sodium azide (see MSDS)* as a preservative.

Precautions

*Due to the sodium azide content a material safety sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling practices.

Product Profile

A minimum dilution of 1:50 is determined by direct immunofluorescent labeling of cultured human fibroblasts or baby hamster kidney (BHK) cells.

In order to obtain best results, it is recommended that each user determine the optimal working dilution for individual applications by titration assay.

F/P Molar Ratio: range is 3-8

Protein Concentration: not less than 1.5 mg/ml by absorbance at 280 nm ($E_{280}^{1\%} = 14.0$)

Storage

For continuous use, store at 2-8 °C for a maximum of one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing is **not** recommended. Storage in "frost-free" freezers is **not** recommended. If slight turbidity occurs upon prolonged

storage, clarify the solution by centrifugation before use.

References

1. Blose, S., et al., J. Cell Biol., **98**, 847 (1984).
2. De la Viña, S., et al., Biochemistry, **27**, 5352 (1988).
3. Masuda, H., et al., J. Cell Biol., **110**, 417 (1990).
4. Schliwa, M., The Cytoskeleton, An Introductory Survey, Springer-Verlag, Wien p. 54 (1985).
5. Mandelkow, E. and Mandelkow, E., In: Guidebook to the Cytoskeletal and Motor Proteins, Kreis, T., and Vale, R. (Eds.), Oxford University Press, Oxford, p 127 (1993).
6. Joshi, H. and Cleveland, D., Cell Mot. Cytosk., **16**, 159 (1990).

Pcs7/00

Sigma brand products are sold through Sigma-Aldrich, Inc.
Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.
Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.
Please see reverse side of the invoice or packing slip.