



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sia.com  
sigma-aldrich.com

## Product Information

**Anti-Human IgG (Fc specific)**  
**Developed in Goat**  
**Delipidized, Whole Antiserum**

Product No. I 8885

### Product Description

Anti-human IgG (Fc specific) is developed in goat using purified Fc fragment of human IgG as the immunogen. The antiserum has been treated to remove lipoproteins.

Specificity for the Fc fragment of human IgG is determined by immunoelectrophoresis (IEP) versus normal human serum and Fc fragment of human IgG. No reactivity is observed versus purified human IgA, IgG (Fab fragment), IgM, Bence Jones kappa, and Bence Jones lambda myeloma proteins.

Identity and purity of the antibody is established by immunoelectrophoresis (IEP). Electrophoresis of the antibody preparation followed by diffusion versus anti-goat IgG results in a single arc of precipitation and against anti-goat serum multiple arcs of precipitation are observed.

### Reagents

Goat Anti-Human IgG is supplied as a liquid containing 0.1% sodium azide as a preservative.

### Precautions

Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

This goat antiserum was maintained at pH 5.0 for 40 minutes to meet U.S.D.A. requirements.

### Product Profile

Each milliliter of antiserum contains 4.0-7.0 mg of specific antibody. Normal human serum is used to determine the antibody concentration by a quantitative precipitation assay.

### Storage

Store at 2-8 °C for up to one month.

For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

PCS 2/00

Sigma brand products are sold through Sigma-Aldrich, Inc.  
Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.  
Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.  
Please see reverse side of the invoice or packing slip.