


SIGMA QUALITY CONTROL TEST PROCEDURE

**Enzymatic Assay of CATALASE
(EC 1.11.1.6)****PRINCIPLE:****CONDITIONS:** T = 25°C, pH = 7.0, A_{240nm}, Light Path = 1 cm**METHOD:** Continuous Spectrophotometric Rate determination**REAGENTS:**

- A. 50 mM Potassium Phosphate Buffer pH 7.0 at 25°C
(Prepare 200 ml in deionized water using Potassium Phosphate, Monobasic, Anhydrous, Sigma Prod. No, P-5379. Adjust to pH 7.0 at 25°C using 1 M KOH.)
- B. 0.036% (w/w) Hydrogen Peroxide Solution (H₂O₂)(Substrate Solution)
(Prepare in Reagent A using Hydrogen Peroxide, 30% (w/w), Sigma Prod. No. H-1009. Determine the A_{240nm} of this solution using Reagent A as a blank. The A_{240nm} should be between 0.550 and 0.520 absorbance units. Add hydrogen peroxide to increase the absorbance and Reagent A to decrease the absorbance.)
- C. Catalase Solution¹
(Immediately before use prepare a solution containing 50 - 100 units per ml in cold Reagent A.)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable cuvettes.

	<u>Test</u>	<u>Blank</u>
Reagent A	-----	3.0
Reagent B (Substrate Solution)	2.9	-----

**Enzymatic Assay of Catalase
(EC 1.11.1.6)**

PROCEDURE: (continued)

Equilibrate to 25°C. Monitor the $A_{240\text{nm}}$ until constant and record the $A_{240\text{nm}}$ using a suitably thermostatted spectrophotometer. Then add:

	<u>Test</u>	<u>Blank</u>
Reagent C (Catalase Solution)	0.1	-----

Immediately mix by inversion and record the time required for the $A_{240\text{nm}}$ to decrease from 0.45 to 0.40 absorbance units.²

CALCULATION:

$$\text{Units/ml enzyme} = \frac{(3.45) (\text{df})}{(\text{min}) (0.1)}$$

3.45 = Corresponds to the decomposition of 3.45 micromoles of hydrogen peroxide in a 3.0 ml reaction mixture producing a decrease in the $A_{240\text{nm}}$ from 0.45 to 0.40 absorbance units

df = Dilution factor

min = Time in minutes required for the $A_{240\text{nm}}$ to decrease from 0.45 to 0.40 absorbance units

0.1 = Volume (in milliliter) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

UNIT DEFINITION:

One unit will decompose 1.0 μmole of H_2O_2 per minute at pH 7.0 at 25°C, while the H_2O_2 concentration falls from

10.3 mM to 9.2 mM. The rate of disappearance of H_2O_2 is followed by observing the rate of decrease in the absorbance at 240 nm.

FINAL ASSAY CONCENTRATION:

In a 3.00 ml reaction mix, the final concentrations are 50 mM potassium phosphate, 0.035% (w/w) hydrogen peroxide, and 5 - 10 units catalase.

**Enzymatic Assay of Catalase
(EC 1.11.1.6)**

REFERENCES:

Beers, R.F. Jr. and Sizer, I.W. (1952) *Journal of Biological Chemistry* **195**, 133-140

Stern, K.G. (1937) *Journal of Biological Chemistry* **121**, 561-572

NOTES:

1. To obtain the proper activity, catalase should be completely dissolved before assaying. Catalase products supplied as a crystalline suspension may require incubation at 37°C upon dilution to obtain complete dissolution.
2. The initial $A_{240\text{nm}}$ will exceed 0.450 absorbance units and start to decrease. Start timing the reaction when the $A_{240\text{nm}}$ reaches 0.450 absorbance units.
3. This assay is based on the cited references.
4. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.