


Anti-p120^{ctn} (P 1870)



Immunoblot

Cytosolic fraction of rat brain extract (Lanes A, C) and ECV304 whole cell extract (Lane B) were separated on SDS-PAGE and blotted with Rabbit Anti-p120^{ctn} (Product No. P 1870). The antibody was developed with Goat Anti-Rabbit IgG-Peroxidase (Product No. A 0545) and chemiluminescent substrate.

Lanes A, B: Antibody dilution 1:2,000

Lane C: Antibody dilution 1:2,000 + 20 μ g/ml p120^{ctn} immunizing peptide