



#501413 Ma Quick Connect
w/#501101 Fe Cap



#500976 Fe Quick Connect
w/#501100 Ma Cap



2.5' #500434 Size 15 C-Flex
Tubing w/Pinch Clamp




2.5' #500434 Size 15 C-Flex
Tubing w/Pinch Clamp



#500138 Spike Port



501796 10L JRH EVA Bag




P1

P2

P3

Double Bagged
w/#501754 JRH Label
& #501403 Irradiation Dot



Each Product Bag Will Be Individually
Double Poly Bagged
& All Tubing Will
Be Coiled



CONFIDENTIAL; THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH
ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

BIOEAZE™ 5-10L EVA ST BAG

Item #1523B