


5' #500145 Size 16 C-Flex Tubing

9" #501603 3/16" x 1/4"
Rigid Polypropylene Tubing

Single Notch Cut

#500123 Tubing Plug

Double Polybagged w/#501754
JRH Label and #501403
Irradiation Dot attached



Each Tubing Set Will Be Individually
Double Poly Bagged
& All Tubing Will
Be Coiled



CONFIDENTIAL; THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH
ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

BIOEAZE™ DISP TUBE OPQ ST

Item #1655R