


5' #501610 1/8" x 1/4" Clear C-Flex Tubing

9" #501603 3/16" x 1/4"  
Rigid Polypropylene Tubing

Single Notch Cut

#500123 Tubing Plug

Double Polybagged w/#501754  
JRH Label and #501403  
Irradiation Dot attached



Each Tubing Set Will Be Individually  
Double Poly Bagged  
& All Tubing Will  
Be Coiled



**CONFIDENTIAL**; THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH  
ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

**BIOEAZE™ DISP TUBE CL ST**

**Item #1656R**