

#501623 Male Luer
w/#501016 Female Cap

Single Poly Bagged
w/#501754 JRH Label &
#501403 Irradiation Dot

Individually Bagged

3' #501610 1/8" x 1/4" Clear
C-Flex Tubing w/Pinch Clamp

X 5

#501413 Ma Quick
Connect w/Vinyl Cap

6" #501947 1/4" x 3/8"
Clear C-Flex Tubing
w/Pinch Clamp

#500138 Spike Port

#501985 2L JRH EVA Bag

P1

P2

P3

Each Product Bag is Individually
Bagged & All Tubing Will
Be Coiled



CONFIDENTIAL: THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

BIOEAZE™ 2L EVA ST BAG ACD 5/PK

Item #1933B