

#501623 Male Luer
w/#501016 Female Cap

Individually Bagged

Single Poly Bagged
w/#501754 JRH Label &
#501403 Irradiation Dot

3' #501610 1/8" x 1/4" Clear
C-Flex Tubing w/Pinch Clamp

#501413 Ma Quick Connect
w/Vinyl Cap

6" #501947 1/4" x 3/8" Clear
C-Flex Tubing w/Pinch Clamp

#500138 Spike Port

#501756 5L JRH EVA Bag

P1 P2 P3

X 5

Each Product Bag is Individually
Bagged & All Tubing Will
Be Coiled



CONFIDENTIAL; THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH
ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

BIOEAZE™ 5L EVA ST BAG ACD 5/PK

Item #1934B