



Single Poly Bagged
w/#501754 JRH Label &
#501403 Irradiation Dot

Individually Bagged



X 5



Each Product Bag is Individually
Bagged & All Tubing Will
Be Coiled

#501413 Ma Quick Connect
w/Vinyl Cap

Same as P1

6" #501947 1/4" x 3/8" Clear
C-Flex Tubing w/Pinch Clamp

#500138 Spike Port

#501726 1L JRH EVA Bag

P1 P2 P3



CONFIDENTIAL; THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH
ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

BIOEAZE™ 1L EVA ST BAG 5/PK

Item #1936B