



Single Poly Bagged w/#501754 JRH Label & #501403 Irradiation Dot

Individually Bagged



X 5



Each Product Bag is Individually Bagged & All Tubing Will Be Coiled


#501413 Ma Quick Connect w/Vinyl Cap

Same as P1

6" #501947 1/4" x 3/8" Clear C-Flex Tubing w/Pinch Clamp

#500138 Spike Port

#501756 5L JRH EVA Bag



**CONFIDENTIAL**: THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

**BIOEAZE™ 5L EVA ST BAG 5/PK**

**Item #1938B**