

#501413 Ma Quick Connect
w/Vinyl Cap


2.5' #501947 1/4" x 3/8" Clear
C-Flex Tubing w/Pinch Clamp

#500123 Tubing Plug


Double Poly Bagged
w/501754 JRH Label &
#501403 Irradiation Dot

3' #501610 1/8" x 1/4" Clear
C-Flex Tubing w/Pinch Clamp

#502102 50L JRH Bulk Bag



Each Item Will Be Individually
Double Poly Bagged & All Tubing Will
Be Coiled



CONFIDENTIAL; THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH
ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

BIOEAZE™ 50L PE ST BAG TER

Item #1946B