

#501091 Ma Quick Connect  
w/#501104 Fe Cap

#501092 Fe Quick Connect  
w/#501103 Ma Cap

Double Poly Bagged  
w/501754 JRH Label &  
#501403 Irradiation Dot


3' #500451 1/2" x 3/4" Silicone  
Tubing w/Pinch Clamp

3' #500451 1/2" x 3/4" Silicone  
Tubing w/Pinch Clamp

#501092 Fe Quick Connect  
w/#501103 Ma Cap

1' #500451 1/2" x 3/4" Silicone  
Tubing w/Pinch Clamp

#501983 JRH 500L Bulk Bag



Each Product Bag Folded, Individually  
Double Bagged & All Tubing Will  
Be Coiled



**CONFIDENTIAL**; THIS IS PROPRIETARY INFORMATION AND SHOULD NOT BE SHARED WITH  
ANY PARTY w/o WRITTEN CONSENT OF JRH BIOSCIENCES.

**BIOEAZE™ 500L PE ST BAG**

**Item #1948B**

P1

P2

P3