



Product Information

Glucose-6-phosphate Dehydrogenase from baker's yeast (*S. cerevisiae*)

Product Number **G 6378**
Storage Temperature -0 °C

Product Description

EC Number: 1.1.1.49
CAS Number: 9001-40-5
Molecular Weight: 128 kDa (gel filtration)¹
Synonym: G-6-P-DH

Glucose 6-phosphate dehydrogenase (G-6-P-DH) is a key regulatory enzyme in the first step of the pentose phosphate pathway. G-6-P-DH oxidizes glucose-6-phosphate in the presence of NADP⁺ to give 6-phosphogluconate. Polyacrylamide gel electrophoresis, activity staining, and anti-yeast G-6-P-DH antibody immunoblotting studies have indicated that G-6-P-DH is a glycoprotein.²

G-6-P-DH has been utilized in assays for nicotinamide adenine dinucleotide⁴ and tissue pyridine nucleotides.⁵ For G-6-P-DH from yeast, the K_m values for glucose 6-phosphate and NADP are 2.0 × 10⁻⁵ M and 2.0 × 10⁻⁶ M, respectively, in Tris buffer (pH 8.0) containing 0.01 M MgCl₂ at 38 °C.⁶

This product is a lyophilized powder that is essentially sulfate-free and contains approximately 20% sodium citrate.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

G-6-P-DH is soluble in water (5 mg/ml), yielding a clear, colorless solution. To ensure maximum stability and recovery of activity, reconstitution with 5 mM citrate, pH 7.5, is recommended. G-6-P-DH can also be dissolved at 1 mg/ml in 5 mM glycine buffer, pH 8, with or without 0.1% BSA. Phosphate buffer will inhibit the enzyme and should not be used.³

Storage/Stability

G-6-P-DH solutions dissolved at 1 mg/ml in 5 mM glycine buffer, pH 8, (with or without 0.1% BSA) or in deionized water at 1 mg/ml can be aliquoted and stored frozen for at least 2 months. However, it is best to subject the aliquots to no more than one freeze/thaw cycle. An aqueous solution can be refrigerated for up to two weeks. For long-term refrigerated storage (> 2 weeks), an ammonium sulfate (3.2 M) suspension of G-6-P-DH is recommended (e.g., Product No. G 7877).

References

1. Andrews, P., The Gel-Filtration Behaviour of Proteins Related to their Molecular Weights over a Wide Range. *Biochem. J.*, **96(3)**, 595-606 (1965).
2. Reilly, K. E., and Allred, J. B., Glucose-6-phosphate Dehydrogenase from *Saccharomyces cerevisiae* Is a Glycoprotein. *Biochem. Biophys. Res. Commun.*, **216(3)**, 993-998 (1995).

3. Domagk, G. F., and Chilla, R.,
Glucose-6-phosphate Dehydrogenase from
Candida utilis. Meth. Enzymol., **XLI-B**, 205-208
(1975).
4. Bernofsky, C., and Swan, M., An Improved Cycling
Assay for Nicotinamide Adenine Dinucleotide.
Anal. Biochem., **53(2)**, 452-458 (1973).
5. Nisselbaum, J. S., and Green, S., A Simple
Ultramicro Method for Determination of Pyridine
Nucleotides in Tissues. Anal. Biochem., **27(2)**,
212-217 (1969).
6. Enzyme Handbook, Vol. I, Barman, T. E.,
Springer-Verlag (New York, 1969), pp. 73-74.

GCY/RXR 11/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.