


### pFLAG-CMV-3 (6.3 kb)



### pFLAG-CMV-4 (6.3 kb)



## Multiple Cloning Site

(pFLAG-CMV-3\* and pFLAG-CMV-4)

FLAG Peptide Sequence														
Met*	Asp	Tyr	Lys	Asp	Asp	Asp	Asp	Lys						
ATG	GAC	TAC	AAA	GAC	GAT	GAC	GAC	AAG	CTT	GCG	GCC	GCG	AAT	TCA
TAC	CTG	ATG	TTT	CTG	CTA	CTG	CTG	TTC	GAA	CGC	CGG	CGC	TTA	AGT
										Hind III				
			Bgl II		EcoR V			Kpn I		Xba I		BamH I		
TCG	ATA	GAT CTG		ATA	TCG	GTA	CCA	GTC	GAC	TCT	AGA	GGA	TCC	CGG
AGC	TAT	CTA	GAC	TAT	AGC	CAT	GGT	CAG	CTG	AGA	TCT	CCT	AGS	GCC

\*For pFLAG-CMV-3, the Met-preprotrypsin leader sequence (PPT LS) precedes the FLAG coding sequence.