



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Monoclonal Anti-Desmosomal Cytokeratin

Mouse Ascites Fluid
Clone DK 80.20

Product No. **C 1041**

Product Description

Monoclonal Anti-Desmosomal Cytokeratin (mouse IgM isotype) is derived from the hybridoma produced by the fusion of mouse myeloma cells and splenocytes from an immunized mouse. Cytoskeleton from a bovine mammary gland epithelial (BMGE) preparation was used as the immunogen. The isotype was determined by double diffusion assay. The product is provided as ascites fluid containing 0.1% sodium azide (see MSDS)* as a preservative.

Monoclonal Anti-Desmosomal Cytokeratin is specific for the cytokeratin peptide 8 by immunoblotting on human and bovine cultured cell lines. Staining shows a typical punctate dot-like pattern consistent with desmosome decoration. The antibody also stains cytokeratin 8 in human tissue extracts and purified cytokeratin polypeptides. The product, when used in immunofluorescent labeling of various bovine and human cultured cells, shows staining of desmosomal patterns along the intercellular contacts of the cells. The product reacts with desmosomes in frozen tissue sections of non-stratified epithelia as follows:

Simple Epithelium

- Small intestine
- Colon
- Pancreas
- Salivary gland (acini, ducts)
- Liver (hepatocytes, bile ducts)

Transitional Epithelium

- Urinary bladder

Pseudostratified Epithelium

- Trachea

The antibody shows no reaction with:

Stratified Epithelium

- Bovine snout
- Tongue
- Cornea
- Esophagus

Non-Epithelia Tissue

- Heart (intercalated discs)

Monoclonal Anti-Desmosomal Cytokeratin may be used for immunohistochemical staining of cultured cells or frozen tissue sections of non-stratified human or bovine epithelia. Aldehyde fixation should be avoided.

Precautions

Due to sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

Product Profile

A minimum working dilution of 1:100 was determined by indirect immunofluorescent staining on frozen sections of bovine tissue.

In order to obtain best results, it is recommended that each individual user determine their working dilution by titration assay.

Storage

For continuous use, store 0-5 °C. For extended storage, solution may be frozen in working aliquots. Repeated freezing and thawing is **not** recommended. If slight turbidity occurs upon prolonged storage, clarify by centrifugation before use.

References

1. Moll, R., et al., Cell, **31**, 11 (1982).
2. Moll, R., et al., Lab. Invest., **54**, 4 (1986).
3. Garrod, D., J. Cell Sci., Suppl. **4**, 221 (1986).

Pcs 6/01

Sigma brand products are sold through Sigma-Aldrich, Inc.
Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications
. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.
Please see reverse side of the invoice or packing slip.