


SIGMA QUALITY CONTROL TEST PROCEDURE

**Enzymatic Assay of
ANGIOTENSIN CONVERTING ENZYME
Sigma Prod. No. A-2580****PRINCIPLE:**

Abbreviations used:

Hip-L-His-L-Leu = Hippuryl-L-Histidyl-L-Leucine

ACE = Angiotensin Converting Enzyme

L-His-L-Leu = L-Histidyl-L-Leucine

CONDITIONS: T = 37 °C, pH = 8.3, A_{214nm}

METHOD: HPLC Stopped Analysis of Products

REAGENTS:

- A. 100 mM Tris HCl Buffer with 300 mM Sodium Chloride, and 10 μM Zinc Chloride, pH 8.3 at 37°C
(Prepare 100 ml in deionized water using Trizma Base, Sigma Prod. No. T-1503, Sodium Chloride, Sigma Prod. No. S-9625, and Zinc Chloride, Sigma Prod. No. Z-4875. Adjust to pH 8.3 at 37°C with 1 M HCl.)
- B. 50 mM Hippuryl-L-Histidyl-L-Leucine Solution (HHL)
(Prepare 1 ml in Reagent A using Hippuryl-His-Leu, Free Base, Sigma Prod. No. H-1635.
PREPARE FRESH.)
- C. 0.02% (w/v) Hippuric Acid Standard Solution (Hipp Std)
(Prepare 1 ml in Reagent A using Hippuric Acid, Free Acid, Sigma Prod. No. H-6375.
Further dilute with Reagent A to standards containing the following concentrations (mg/ml):
0.05, 0.075, 0.10, and 0.15.)
- D. Acetonitrile
(Use Acetonitrile, Sigma Stock No. 27,071-7.)
- E. 0.08% (v/v) Phosphoric Acid, pH 2.5
(Prepare 100 ml in deionized water using Phosphoric Acid, Sigma Prod. No. P6560. The pH of a 0.08% (v/v) is approximately 2.5.)

**Enzymatic Assay of
ANGIOTENSIN CONVERTING ENZYME
Sigma Prod. No. A-2580**

REAGENTS: (continued)

- F. Angiotensin Converting Enzyme Solution
(Immediately before use, prepare a solution containing 0.05-0.10 unit/ml of Angiotensin Converting Enzyme in cold Reagent A.)

PROCEDURE:

Step 1:

Pipette (in milliliters) the following reagents into suitable microcentrifuge tubes:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	0.08	0.09
Reagent F (Enz Sol)	0.01	-----

Equilibrate to 37°C. Then add:

Reagent B (HHL)	0.01	0.01
-----------------	------	------

Incubate for exactly 30 minutes at 37°C.

Terminate the reaction by heating at 100°C for 4 minutes and then microcentrifuge at 15,000 g for 10 minutes. Transfer 0.08 ml of both the Test and Blank to HPLC vials.

Step II:

HPLC Analysis of Products

The substrate and products are separated by reverse phase HPLC chromatography as follows:

1. Column: Supelcosil LC-18, Supelco Prod. No. 5-8298,
Particle size: 5 µm, 25 cm x 4.6 mm.
2.

Mobile Phase	Time (min)		
Reagent D (Acetonitrile)	0	15	25
Reagent E (Phosphoric Acid)	4.5	30	30
	95.5	70	70
3. Pressure: 2150 PSI², Flow rate 1.5 ml/min, Detection: 214nm, Sample Injection volume: 20 µl
4. Inject Blank and then Standards of Reagent C (Hipp Std) and samples. A comparison can then be made between the standard curve of hippuric acid and hippuric acid generated from the sample reaction. Convert mg of hippuric acid to µmoles.

**Enzymatic Assay of
ANGIOTENSIN CONVERTING ENZYME
Sigma Prod. No. A-2580**

CALCULATIONS:

$$\text{Units/ml enzyme} = \frac{(\mu\text{moles of hippuric acid})(df)(0.1)}{(0.01)(30)}$$

0.1 = Total volume (in milliliter) of enzyme assay
0.01 = Volume (in milliliter) of enzyme used
30 = Time (in minutes of assay)

$$\text{Units/mg solid} = \frac{\text{Units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{Units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

FINAL ASSAY CONDITIONS:

In a 0.10 ml reaction mix, the final concentrations are 100 mM Tris, 300 mM sodium chloride, 10 μ M zinc chloride, 5 mM Hippuryl-L-histidyl-L-leucine and 0.0005-0.001 unit angiotensin converting enzyme.

UNIT DEFINITION:

One unit will produce 1.0 μ mole of hippuric acid from hippuryl-His-Leu per min in 100 mM Tris HCl, 300 mM NaCl, and 10 μ M ZnCl₂ at pH 8.3 at 37°C.

REFERENCE:

Hooper, N.M., and Turner, A.J. (1987) Biochemical Journal **241**, 625-633

NOTES:

1. This assay is not to be used to assay Angiotensin Converting Enzyme, Sigma Prod. No. A-6778.
2. This pressure may vary.
3. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.
4. This assay is based on the cited reference.

Sigma warrants that the above procedure information is currently utilized at Sigma and that all Sigma-Aldrich, Inc. products conform to the information in this and other Sigma-Aldrich, Inc. publications. Purchaser must determine the suitability of the information and product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.